macro

Gebrauchsanleitung • Operating Manual • Mode d'emploi
Instrucciones de manejo • Istruzioni per l'uso

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Tel.: +1-860-767 2562
Fax: +1-860-767 2563
www.brandtech.com

Indien:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park, Powai
Mumbai - 400 076 (India)

Tel.: +91 22 42957790
Fax: +91 22 42957791
E-Mail: info@brand.co.in
www.brand.co.in

China:

BRAND (Shanghai) Trading Co., Ltd.
Guangqi Culture Plaza
Room 506, Building B
No. 2899, Xietu Road
Shanghai 200030 (P.R. China)

Tel.: +86 21 6422 2318
Fax: +86 21 6422 2268
E-Mail: info@brand.cn.com
www.brand.cn.com

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	Page
Safety Instructions	20
Function and Limitations of Use	21
Operating Elements	22
Pipetting	23
Cleaning · Autoclaving	26
Replacing the Filter and Cleaning the Adapter	26
Replacing the Valve System	28
Leak Test	29
Ordering Data · Spare Parts	30
Troubleshooting	31
Repairs	32
Markings on the product	32
Contact addresses	33
Warranty Information	34

Safety Instructions

This instrument may sometimes be used with hazardous materials, operations, and equipment. It is beyond the scope of this manual to address all of the potential safety risks associated with its use in such applications. It is the responsibility of the user of this instrument to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

Please read the following carefully!

1. Every user must read and understand this operating manual before operation.
2. Follow general instructions for hazard prevention and safety instructions; e.g., wear protective clothing, eye protection and gloves.
3. Observe all specifications provided by reagent manufacturers.
4. Use the instrument only for pipetting liquids, with strict regard to the defined limitations of use. Observe operating exclusions (s. page 21)! If in doubt, contact the manufacturer or supplier.
5. Always use the instrument in such a way that neither the user nor any other person is endangered. Avoid splashes. Only use suitable vessels.
6. Never use force on the instrument!
7. Use only original manufacturer's accessories and spare parts. Do not attempt to make any technical alterations. Do not dismantle the instrument any further than is described in the operating manual!
8. Always check the instrument for visible damage before use. If there is a sign of a potential malfunction, immediately stop pipetting. Consult the 'Troubleshooting' section of this manual (see page 31), and contact the manufacturer if needed.

The instrument is designed to assist the filling and dispensing of graduated and volumetric pipettes of glass or plastic in the volume range of 0.1 ml to 200 ml with a suction tube outer diameter < 9.2 mm for measuring liquids. If the instrument is used correctly, the pipetted liquid will only contact the pipette.

Limitations of Use

The instrument is designed for pipetting liquids, observing the following physical limits:

- from $+10\text{ }^{\circ}\text{C}$ to $+40\text{ }^{\circ}\text{C}$ (from $50\text{ }^{\circ}\text{F}$ to $104\text{ }^{\circ}\text{F}$)
(instrument and reagent)
- Vapor pressure up to max. 500 mbar.
Aspirate slowly above 300 mbar, in order to prevent the liquid from boiling.
- Density up to 2.4 g/cm^3

Operating Exclusions

Never use the instrument with liquids whose vapors attack the materials silicone or PP.


The instrument is not designed for Pasteur pipettes.

Storage Conditions

Store the instrument and accessories only in cleaned condition in a cool and dry place.

Storage temperature: $-20\text{ }^{\circ}\text{C}$ to $+50\text{ }^{\circ}\text{C}$ ($-4\text{ }^{\circ}\text{F}$ to $+122\text{ }^{\circ}\text{F}$).

Operating Elements

1. Compress the suction bellows

Before attaching the pipette, squeeze the suction bellows.

2. Attach the pipette

Hold the pipette as near to its upper end as possible, and carefully insert it into the adapter until it fits tightly.

Once the pipette has been securely attached, always hold the instrument in a vertical position, tip down.

Warning!


Be sure that the pipette fits tightly into the adapter. Never use force. Thin pipets are particularly liable to break. Avoid the risk of injury!

3. Fill the pipette

Immerse the pipette tip into the liquid. Press the pipetting lever slowly upwards. Fill the pipette so that the level of liquid is slightly above the required mark.

Attention!

Please take care that no liquid enters into the instrument. This would impair the filtering function and reduce the suction capacity. If liquid does enter the instrument, exchange the filter (see pages 26 - 27).

Note:

The further up the lever is pressed, the stronger the suction becomes and the quicker the meniscus rises.

In the case of large pipettes > 50 ml the vacuum contained in the suction bellows is not sufficient to draw in all the liquid at once. Therefore, squeeze the suction bellows again and continue drawing up liquid.

4. Adjust the volume

Use suitable lint-free tissue to wipe the pipette tip, if necessary. Press the pipetting lever down slowly, until precise adjustment of the meniscus has been achieved.

5. Dispense the pipette

Hold the collecting vessel in an inclined position. Place the pipette tip against the inner vessel wall and press the pipetting lever down to dispense.

5.1 In the case of pipettes with a waiting time

(Imprint e.g., 'Ex + 5 s'):

- As soon as the meniscus in the pipette comes to a standstill, start waiting time as is indicated on the pipette (here: 5 s).
- Wipe the pipette tip a few millimeters upward along the wall of the vessel.

5.2 In the case of blow-out pipettes

(Imprint 'ausblasen – blow out'):

- As soon as the meniscus in the pipette tip comes to a standstill, press the small blow-out bellows once.
- Wipe the pipette tip a few millimeters upward along the wall of the vessel.

6. After pipetting

Hold the pipette as its extreme upper end, and gently twist and pull it out of the adapter.

When properly used, the instrument is maintenance-free. Before autoclaving, the instrument must be cleaned carefully. Wipe the housing with a damp cloth. The instrument is autoclavable at 121 °C (250 °F), 2 bar absolute (30 psi) with a holding time of at least 15 minutes according to DIN EN 285.

Replacing the Filter and Cleaning the Adapter

Only necessary if a pipette has been overfilled or the suction capacity is impaired.

Warning!

Wear protective gloves and eye protection.
Avoid the risk of contamination!

1. Pull out the pipette.
2. Unscrew the adapter support.
3. Pull out the membrane filter.
4. Remove the adapter from the adapter support by pushing it upwards.
5. To clean adapter, carefully rinse it by using an appropriate compatible solvent (e.g., water) in a wash bottle. Blow out the liquid and allow to dry completely.

6. Place the new membrane filter with the thicker end facing downwards into the adapter (with the smaller conical side up).

7. Screw in the adapter support hand tight.

8. The entire instrument incl. membrane filter is autoclavable at 121 °C (2 bar) according to DIN EN 285.

For autoclaving remove suction bellow.

Note:

Not every membrane filter is suitable. Only use genuine manufacturer's recommended accessories. Filters of a pore size of 3 µm are a standard equipment of the instrument. The membrane filter is autoclavable up to 5 times at 121 °C. It is the user's responsibility to ensure effective autoclaving.

Replace the Valve System

1. Unscrew the suction bellows and the adapter support.

2. Take off the upper part of the housing. Carefully pull out the valve system and replace it.

3. Reassemble the instrument in reverse order and carry out leak test.

Leak Test

When the meniscus has been set, liquid should not drip out of the pipette until the valves are activated.

Should the pipette drip, see chapter 'Troubleshooting' (see page 31).

macro pipette controller,

with spare membrane filter 3 µm and operating manual.

Color*	Cat. No.
grey	262 00
green	262 01
blue	262 02
magenta	262 03

* Colors refer to accent colors on instrument adapter support and locking ring.

Adapter support,

PP, length 49 mm

Color	Cat. No.
grey	262 20
green	262 23
blue	262 21
magenta	262 22

Spare membrane filter 3 µm,

non-sterile, 1 pc. in blister pack

Cat. No.	260 52
----------	--------

Spare membrane filter 3 µm,

non-sterile, 10 pcs., in PE-bag

Cat. No.	260 56
----------	--------

Suction bellows with locking ring,

Color	Cat. No.
grey	262 25
green	262 28
blue	262 26
magenta	262 27

Silicone adapter,

length 44 mm

Cat. No.	261 46
----------	--------

Valve system

Cat. No.	261 28
----------	--------

Trouble	Possible Cause	Action to be taken
Suction capacity impaired	Filter dirty or valve system damaged	Replace filter, see pages 26-27
Pipette drips	Filter not properly positioned or valve system damaged	Insert filter properly or replace valve system, see page 28
Pipette not held properly	Silicone adapter damaged	Unscrew the adapter support of the instrument, replace silicone adapter
Compressed suction bellows take in air automatically	Suction bellows or valve system damaged	Replace the locking ring with the suction bellows or the valve system

Note:

Often visual control of the instrument is sufficient to discover a faulty device which can easily be replaced by a spare part. The advantages of do-it-yourself repair are obvious: It saves time and money.

Markings on the product

Markings	Meaning
	Follow the instructions listed on the device, the accessory parts and in the operating manual.
	Observe the operating manual.

Return for repair

Important!

Transporting of hazardous materials without a permit is a violation of federal law.

- Clean and decontaminate the instrument carefully.
- It is essential always to include an exact description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

Outside the U.S. and Canada:

- Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de.

In the U.S. and Canada:

- Contact BrandTech Scientific, Inc. and obtain authorization for the return **before** sending your instrument for service.
- Return only cleaned and decontaminated instruments, with the Return Authorization Number prominently displayed on the outside of the package to the address provided with the Return Authorization Number.

BRAND GMBH + CO KG

Otto-Schott-Straße 25
97877 Wertheim (Germany)

Tel.: +49 9342 808-0
Fax: +49 9342 808-98000
E-Mail: info@brand.de
www.brand.de

USA and Canada:

BrandTech® Scientific, Inc.
11 Bokum Road
Essex, CT 06426-1506 (USA)
Tel.: +1-860-767 2562
Fax: +1-860-767 2563
www.brandtech.com

India:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park, Powai
Mumbai - 400 076 (India)
Tel.: +91 22 42957790
Fax: +91 22 42957791
E-Mail: info@brand.co.in
www.brand.co.in

China:

BRAND (Shanghai) Trading Co., Ltd.
Guangqi Culture Plaza
Room 506, Building B
No. 2899, Xietu Road
Shanghai 200030 (P.R. China)
Tel.: +86 21 6422 2318
Fax: +86 21 6422 2268
E-Mail: info@brand.cn.com
www.brand.cn.com

Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the instrument or the consequences of normal wear and tear especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow the instructions of the operating manual. We are not liable for damage resulting from any actions not described in the operating manual or if non-original spare parts or components have been used.

U.S. and Canada:

Information for warranty please see www.brandtech.com.

Subject to technical modification without notice. Errors excepted.

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