HandyStep® S

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Operating Manual · 23

Mode d'emploi · 43

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	Page
Safety Instructions	25
Function and Limitations of Use	26
Operating and Control Elements	27
Inserting the PD-Tip	28
Setting the Volume	29
Filling the PD-Tip	30
Repetitive Dispensing	31
Ejecting the PD-Tip	32
Checking the Volume	33
Accuracy Table	35
Servicing · Cleaning · UV sterilization	36
Ordering Information · Accessories	37
Troubleshooting	38
Repairs	39
Contact addresses	40
Calibration Service	41
Warranty Information	42
Disposal	42

Please read the following carefully!

This instrument may sometimes be used with hazardous materials, operations, and equipment. It is beyond the scope of this manual to address all of the potential safety risks associated with its use in such applications. It is the responsibility of the user of this pipette to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

- 1.** Every user must read and understand this operating manual prior to using the instrument and observe these instructions during use.
- 2.** Follow general instructions for hazard prevention and safety instructions; e.g., wear protective clothing, eye protection and gloves. When working with infectious or other hazardous samples, all appropriate regulations and precautions must be followed.
- 3.** Observe the reagent manufacturers' information.
- 4.** Only use the instrument for dispensing liquids that conform to the specifications defined in the limitations of use and operating limitations. Observe operating exclusions (see page 26). If in doubt, contact the manufacturer or supplier.
- 5.** Always use the instrument in such a way that neither the user nor any other person is endangered. Avoid splashes. Use only suitable vessels.
- 6.** Avoid touching the tip orifices when working with hazardous samples.
- 7.** Never use force on the instrument!
- 8.** Do not attempt to make any technical alterations. Do not dismantle the instrument.
- 9.** Before use, check the instrument for visible damages. If there is a sign of a potential malfunction, immediately stop dispensing. Consult the 'Troubleshooting' section of this manual (see page 38), and contact the manufacturer if needed.

Function and Limitations of Use

The HandyStep® S is a repetitive pipette for quick and simple repetitive dispensing of liquids. Volumes from 2 µl to 5 ml can be dispensed with the highest precision and accuracy through combinations of the PD-Tip size used and adjustment of the stroke setting on the stroke-setting wheel.

Limitations of Use

This instrument is designed for dispensing samples, observing the following limits:

- use between +15 °C and +40 °C (59 °F and 104 °F)
(of instrument and reagents – other temperatures on request)
- vapor pressure up to 500 mbar
- viscosity: 20 mPa s with 50 ml PD-Tips
260 mPa s with 5 ml PD-Tips
977 mPa s with 1 ml PD-Tips

Operating Exclusions

During proper operation, the dispensed sample only comes in contact with the tip and not with the HandyStep® S.


The user has to ensure the compatibility of the instrument with the intended application. In case of third-party dispenser tips, the user must check their compatibility and perform a function and volume test before use.

The instrument should not be used for dispensing liquids that attack polypropylene (e.g., PD-Tip cylinder), polyethylene (PD-Tip piston), LCP (PD-Tip piston 0.1 ml), or PC/PBT and polycarbonate.

Note:

Instrument and tips are not autoclavable (sterile PD-Tips, see page 37). PD-Tips are disposables. For disinfection of the instrument you may use commercially available disinfecting solutions as long as they do not attack PC/PBT and PC.

Operating and Control Elements

The volume table on the instrument backside is printed on both sides. To change the table, simply push the clip up, turn the table around, and push the clip back down again.

Front side:
stroke settings

Backside:
volume table

Inserting the PD-Tip

- Use of the 25 ml and 50 ml PD-Tips requires the reusable adapter supplied with the PD-Tips. This is coupled to the PD-Tips with a bayonet lock, and can be detached after use.
- PD-Tips are disposables.

- 1.** Push the locking/filling lever to the lower stop.

- 2.** Swing out the locking/filling lever.

- 3.** Insert PD-Tip straight and upright from below.

- 4.** Swing back the locking/filling lever.

Adapter

The possible dispensing volumes of the HandyStep® S are listed in tables. These are found beneath the transparent clip on the backside of the instrument. There are 9 half-integral stroke setting positions, from 1 to 5.

1. Look for the desired dispensing volume (Volume) in the table.
2. Use the stroke setting wheel to set to the proper stroke number (Setting) for the volume that corresponds to the PD-Tip (Tip size).
3. Dispensing steps: depending on the PD-Tip capacity (ml), the same desired dispensing volume per step can be delivered a different number of times (see example below). The smaller the number of steps, the higher the accuracy of the dispensed volume.

Example:

Desired dispensing volume per step: **200 µl**
 PD-Tip sizes to use:

2.5 ml = stroke setting **4** = max. **11** dispensing steps

5 ml = stroke setting **2** = max. **24** dispensing steps

10 ml = stroke setting **1** = max. **49** dispensing steps

Setting	Tip size (ml)										
	0.1	0.5	1	1.25	2.5	5	10	12.5	25	50	
1	2	10	20	25	50	100	200	250	500	1000	49
1.5	3	15	30	37.5	75	150	300	375	750	1500	32
2	4	20	40	50	100	200	400	500	1000	2000	24
2.5	5	25	50	62.5	125	250	500	625	1250	2500	19
3	6	30	60	75	150	300	600	750	1500	3000	15
3.5	7	35	70	87.5	175	350	700	875	1750	3500	13
4	8	40	80	100	200	400	800	1000	2000	4000	11
4.5	9	45	90	112.5	225	450	900	1125	2250	4500	10
5	10	50	100	125	250	500	1000	1250	2500	5000	9
Volume (µl)											

Important!

The volumes listed in the table are the only options available. Since the first dispensing step must be discarded, the number given in the table is always 1 step less.

Filling the PD-Tip

Aspirating Liquid


Note:

We recommend that a new PD tip has to be primed before use. For this, fill the tip at a minimum with the liquid and then empty it again. Small air bubbles in the area of the piston after priming do not affect the results.

1. Push the locking/filling lever to the lower stop.
2. Immerse the orifice of the PD-Tip vertically into the liquid to a depth of 3 - 10 mm.
3. Raise the locking/filling lever slowly to the upper stop to avoid the formation of air bubbles.


Note:

Any small air bubble near the piston can be ignored since a locking mechanism prevents the residual volume from being dispensed inadvertently after the final complete dispensing step.

Dispensing Operation

1. Check the volume setting once again. (Compare with the settings in the volume table.)
2. Wipe off any liquid adhering to the outside of the PD-Tip orifice with a lint-free tissue.
3. Place the orifice of the PD-Tip against the wall of the vessel.
4. Dispense the liquid by pressing down the dispensing lever completely, then let it slide all the way back up.
5. Take care to dispense smoothly and evenly.

Important!

The first dispensing step must be discarded!

Ejecting the PD-Tip

Operation

1. Hold the HandyStep® S over a vessel.
2. Empty the PD-Tip by pushing the locking/filling lever to its lower stop.
3. Swing out the locking/filling lever (the piston will then be decoupled).
4. Hold the HandyStep® S over the waste container. Press the ejection key down, releasing the cylinder. The PD-Tip will be ejected.


Warning!

Residues of the medium might be present on the PD-Tip orifice. The PD-Tip should only be ejected in a manner that does not pose a hazard to the user or other persons.

Note:

The tightness of the seal of the PD-Tips is not guaranteed for repeated dispensing of high viscosity media.

Depending on use, we recommend that gravimetric testing of the instrument be carried out every 3-12 months. This time frame should be adjusted to individual requirements. Gravimetric volume testing according to DIN EN ISO 8655-5 is performed as follows:

1. Preparation of the instrument

Fitting tip. This test can be carried out with PD-Tips of any size. However, size 5 ml is most commonly used for this purpose.

2. Carry out the test

- a) Adjust the stroke setting of HandyStep® S to Step 5 (10 % of nominal volume).
- b) Filling the PD-Tip, immersing it vertically into the testing liquid.
- c) Discard the first step; it only serves to align the internal mechanism.
- d) Place the weighing vessel (containing a small amount of deionized water) on a balance and tare the balance.
- e) Dispense the second step into the weighing vessel. For this, push down the dispensing lever at a steady rate to the stop and hold it there. Then wipe the tip on the vessel wall over a distance of approx. 10 mm.
- f) Enter weighed value into test record.
- g) Repeat points D to F 10 times.
(For the 10th volume test, the 5 ml PD-Tip with a stroke setting of 5 has to be filled a second time – steps A through C.)
- h) Repeat the same testing procedure at stroke settings 3 (6 % of nominal volume) and 1 (2 % of nominal volume).
- i) This results in a total of the 30 weighed values needed.

Checking the Volume

Calculation (for nominal volume)

x_i = Weighing results

n = Number of weighings

Z = Correction factor

(e.g., 1.0029 $\mu\text{l}/\text{mg}$
at 20 °C, 1013 hPa)

$$\text{Mean value } \bar{x} = \frac{\sum x_i}{n}$$

$$\text{Mean volume } \bar{V} = \bar{x} \cdot Z$$

Accuracy*

$$\mathbf{A\%} = \frac{\bar{V} - V_0}{V_0} \cdot 100$$

V_0 = Nominal volume

Coefficient of Variation*

$$\mathbf{CV\%} = \frac{100 \cdot s}{\bar{V}}$$

Standard Deviation

$$\mathbf{s} = Z \cdot \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

*) = Calculation of accuracy (A%) and variation coefficient (CV%):
A% and CV% are calculated according to the formulas for statistical control.

Note:

Testing instructions (SOPs) and a demo version of the EASYCAL™ 4.0 calibration software are available for download at www.brand.de.

Technical Data

PD-Tips from BRAND, 20 °C 'Ex'

PD-Tip size	Volume range (µl)	A* ≤ ± % Stroke setting ≅ % of nominal volume			CV* ≤ % Stroke setting ≅ % of nominal volume		
		1 ≅ 2%	3 ≅ 6%	5 ≅ 10%	1 ≅ 2%	3 ≅ 6%	5 ≅ 10%
		0.1 ml	2 - 10	4.0	2.4	1.6	6.0
0.5 ml	10 - 50	2.5	1.5	1.0	2.5	1.5	1.0
1.0 ml	20 - 100	2.5	1.5	1.0	2.0	1.2	0.8
1.25 ml	25 - 125	2.5	1.4	0.9	2.0	1.1	0.7
2.5 ml	50 - 250	1.8	1.1	0.7	1.5	0.9	0.6
5.0 ml	100 - 500	1.8	1.1	0.7	1.5	0.9	0.7
10.0 ml	200 - 1000	1.8	1.1	0.7	2.0	1.2	0.8
12.5 ml	250 - 1250	1.8	1.1	0.8	3.2	2.0	1.4
25.0 ml	500 - 2500	1.5	0.9	0.6	3.0	1.5	1.0
50.0 ml	1000 - 5000	1.5	0.8	0.5	5.0	1.8	1.2

A* = Accuracy, CV* = Coefficient of variation

Error limits refer to the dispensing volume set relative to the PD-Tip, obtained when instrument, PD-Tip and distilled water are equilibrated at ambient temperature (20 °C/68 °F) and with smooth operation. The testing is according to DIN EN ISO 8655-5. The nominal volume is the maximum volume printed on the PD-Tip.

Note:

The device is marked in accordance with the German Measurement and Calibration Act as well as the Measurement and Calibration Regulation:

DE-M 18

Character string: DE-M (DE for Germany), framed by a rectangle, as well as the last two digits of the year in which the marking was affixed (here: 2018).

Important!

The HandyStep® S can also be operated with compatible tips from other manufacturers. The user must check their compatibility and perform a function and volume test before use!

Serviceing · Cleaning · UV sterilization

HandyStep® S is factory calibrated and maintenance-free.
The instrument must not be dismantled!

In cases of external contamination, we recommend cleaning
the instrument using a wipe with water or isopropanol.

UV sterilization

The unit can withstand the usual
output of a UV sterilization lamp. The
effects of the UV may cause some
color change.

Important!

Avoid liquid getting inside the instrument.
In case this happens, let the instrument
be opened and cleaned by an autho-
rized service center, or send it back to
BRAND.

The rack mount is provided with
an adhesive strip for mounting. It
can be used without the strip in the
Transferpette® S 6-position bench-
top rack.

HandyStep® S

pack of 1, 1 shelf/rack mount,
3 PD-Tips (0.1 ml, 1 ml, 10 ml).

Cat. No. 7051 10

Shelf/rack mount
pack of 1.

Cat. No. 7051 30

PD-Tips (Positive Displacement Tips)

Capacity	Pack of	non sterile Cat. No.	sterile** Cat. No.	BIO-CERT® Cat. No.
0.1 ml	100	7024 02	7024 04	7026 83
0.5 ml	100	7023 70	7023 84	7026 84
1.0 ml	100	7024 06	7024 36	7026 85
1.25 ml	100	7023 72	7023 86	7026 86
2.5 ml	100	7023 74	7023 88	7026 88
5.0 ml	100	7023 76	7023 90	7026 90
10.0 ml	100	7024 07	7024 38	7026 91
12.5 ml	100	7023 78	7023 92	7026 92
25.0 ml*	50/25***	7023 80	7023 94	7026 94
50.0 ml*	25	7023 82	7023 96	7026 96

* incl. 1 adapter, ** sterile/free of endotoxin – single wrapped,

*** PD-Tips 25 ml: non sterile pack of 50 / sterile pack of 25

Adapter for 25 and 50 ml

PD-Tips, PP, autoclavable

non sterile, pack of 10

Cat. No. 7023 98

BIO-CERT®, pack of 5

Cat. No. 7023 99

PD-Tip Set

20 PD-Tips each of 0.5, 1.0, 1.25,
2.5, 5.0, 10 and 12.5 ml size

Cat. No. 7023 68

PD-Tips BIO-CERT® are sterile, free of endo-
toxins, DNA, RNase and ATP.

Troubleshooting

Problem	Possible cause	Corrective action
PD-Tip piston cannot be locked after the PD-Tip has been inserted.	Locking/filling lever is not pushed down completely and is not swung out.	First press the ejection key and remove the PD-Tip, then push the locking/filling lever all the way down and swing it out completely.
Undefined dispensing volume	The stroke setting wheel is not properly engaged.	Ensure that the stroke setting wheel is securely locked into place.
Air bubble beneath the PD-Tip piston becomes larger.	Leaking PD-Tip	Replace PD-Tip
A and/or CV exceed tolerances.	Leaking PD-Tip	Replace PD-Tip
	The dispenser tips are difficult to move.	Replace the PD-Tip or use other more pliant dispenser tips.
	The dispensing lever has not been pressed completely and evenly.	Press the dispensing lever evenly to stop position.
	The first dispensing step was not discarded.	Discard the first dispensing step.
The PD-Tip piston cannot be decoupled after emptying completely.	The PD-Tip has been inserted improperly.	Press the ejection key, push down the locking/filling lever completely, and swing it out.

Return for repair

Important! Transporting of hazardous materials without a permit is a violation of federal law.

- Clean and decontaminate the instrument carefully.
- It is essential always to include an exact description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

Outside the U.S. and Canada:

- Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de.

In the U.S. and Canada:

- Contact BrandTech Scientific, Inc. and obtain authorization for the return **before** sending your instrument for service.
- Return only cleaned and decontaminated instruments, with the Return Authorization Number prominently displayed on the outside of the package to the address provided with the Return Authorization Number.

BRAND GMBH + CO KG

Otto-Schott-Straße 25
97877 Wertheim (Germany)

Tel.: +49 9342 808-0

Fax: +49 9342 808-98000

E-Mail: info@brand.de

www.brand.de

USA and Canada:

BrandTech® Scientific, Inc.
11 Bokum Road
Essex, CT 06426-1506 (USA)

Tel.: +1-860-767 2562

Fax: +1-860-767 2563

www.brandtech.com

India:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park, Powai
Mumbai - 400 076 (India)

Tel.: +91 22 42957790

Fax: +91 22 42957791

E-Mail: info@brand.co.in

www.brand.co.in

China:

BRAND (Shanghai) Trading Co., Ltd.
Guangqi Culture Plaza
Room 506, Building B
No. 2899, Xietu Road
Shanghai 200030 (P.R. China)

Tel.: +86 21 6422 2318

Fax: +86 21 6422 2268

E-Mail: info@brand.cn.com

www.brand.cn.com

ISO 9001 and GLP-guidelines require regular examinations of your volumetric instruments. We recommend checking the volume every 3-12 months. The interval depends on the specific requirements on the instrument. For instruments frequently used or in use with aggressive media, the interval should be shorter. The detailed testing instruction can be downloaded on www.brand.de or www.brandtech.com.

BRAND also offers you the possibility to have your instruments calibrated by the BRAND Calibration Service or the BRAND-owned DAkkS Calibration Service. Just send in the instruments to be calibrated, accompanied by an indication of which kind of calibration you wish. Your instruments will be returned within a few days together with a test report (BRAND Calibration Service) or with a DAkkS Calibration Certificate. For further information, please contact your dealer or BRAND. Complete ordering information is available for download at www.brand.de (see Technical Documentation).

Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the instrument or the consequences of normal wear and tear especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow the instructions of the operating manual. We are not liable for damage resulting from any actions not described in the operating manual or if non-original spare parts or components have been used.

U.S. and Canada:

Information for warranty please see www.brandtech.com.

Disposal

For the disposal of instruments and PD-Tips, please observe the relevant national disposal regulations.

Subject to technical modification without notice. Errors excepted.

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