[Gebrauchsanleitung](#) | [Operating manual](#) | [Mode d'emploi](#) |
[Instrucciones de manejo](#) | [Istruzione](#) | [Instruções de utilização](#) | [操作手册](#) |
[Руководство по эксплуатации](#) | [사용 지침](#) | [Bruksanvisning](#) | [使用説明書](#) |
[Használati utasítás](#) | [Návod k použití](#) | [Gebruiksaanwijzing](#) |
[Instrukcja użytkowania](#) | [Kullanım Talimatları](#)

Transferpette® S

Pipettierhelfer | Pipette controllers

Impressum

BRAND GMBH + CO KG

Otto-Schott-Str. 25

97877 Wertheim (Germany)

T +49 9342 808 0

F +49 9342 808 98000

info@brand.de

www.brand.de

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Table of contents

1	Introduction	34
1.1	Scope of delivery.....	34
1.2	Use the instructions for use.....	34
2	Safety regulations.....	35
2.1	Safety regulations	35
2.2	Purpose	35
2.3	Limitations of use.....	35
2.4	Application restrictions	36
2.5	Operating exclusions	36
3	Functional and operational components	37
4	Pipetting	38
5	Checking the volume	41
6	Accuracy table	42
7	Adjustment – Easy Calibration.....	44
8	Disinfection/autoclaving.....	46
8.1	Autoclaving.....	46
8.2	UV sterilization	46
8.3	PE filter	46
9	Maintenance.....	47
9.1	Check pipetting unit	47
9.2	Disassembling/cleaning (up to 1,000 µl)	48
9.3	Disassembly/cleaning (2–10 ml)	49
10	Troubleshooting	50
11	Product markings	51
12	Order Information.....	52
12.1	Order info/accessories.....	52
12.2	Spares.....	53
12.3	Additional accessories.....	55
13	Repairs.....	56
13.1	Sending for repair	56
14	Calibration service	58
15	Information about your laboratory in- strument	59
16	Warranty	60
17	Disposal	61

1 Introduction

1.1 Scope of delivery

Transferpette® S, adjustable/fixed-volume, DE-M marking, supplied with quality certificate, shelf mount, and silicone grease.

1.2 Use the instructions for use

- Carefully read the operating manual before using the device for the first time.
- Keep the user manual easily accessible. It is part of the device.
- Be sure to include the operating manual if you transfer possession of this device to a third party.
- You can find up-to-date versions of the operating manual on our website: www.brand.de.

1.2.1 Signal words and their meaning

Signal word	Meaning
⚠ WARNING or ⚠ WARNING! ...	WARNING indicates a dangerous situation that, if not avoided, could result in death or serious injury.
⚠ CAUTION or ⚠ CAUTION! ...	CAUTION indicates a hazardous situation that, if not avoided, could result in moderate or minor injury.
NOTICE or NOTICE ...	NOTE is used to address actions that are not related to physical injury. Example: Possible property damage.

1.2.2 Presentation of descriptions of actions

Format	Meaning
1. Task	Indicates a task.
a., b., c.	Indicates the individual steps of a task.
>	Indicates a prerequisite for a task.
⇒	Indicates a result of a completed task.

2 Safety regulations

2.1 Safety regulations

Please read carefully!

The laboratory device Transferpette® S can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

- Every user must read and observe the user manual enclosed with this laboratory device before and during use. The device may only be used by trained and qualified personnel.
- Follow the general hazard instructions and safety regulations, e.g. wear protective clothing, eye protection and protective gloves.
- When working with infectious or hazardous samples/media (e.g. hazardous substances), the general safety rules in the laboratory must be observed and regulations for handling the samples/media must be followed. The information provided by the media manufacturer (e.g. safety data sheets) must be observed.
- Use the laboratory device only for pipetting or dispensing media within the defined limitations and restrictions of use. Comply with the operating exclusions.
- If working with inflammable media, make sure to take the necessary precautions to avoid the buildup of static electricity (e.g., do not dispense into plastic vessels and do not wipe instruments with a dry cloth). Do not use the device in potentially explosive atmospheres. If in doubt, contact the manufacturer or supplier.
- Always check that the instrument is in proper working condition before use. If device malfunctions are indicated (e.g. sluggish pistons, leaks or power supply faults), stop working with the device immediately and refer to the "Troubleshooting" section in the user manual. Contact the manufacturer, if necessary.
- Always perform work in a manner that does not endanger users or other people. Avoid splashes. Only use suitable vessels. Never use unnecessary force when operating, cleaning or servicing the laboratory device.
- If the laboratory device is supplied with power through a power adapter, batteries or accumulators, the components and connections on the device must be regularly checked for proper working condition. Do not operate the laboratory device and its accessories in an unprotected, humid or wet environment.
- Do not make any technical modifications. Use only original replacement parts of the manufacturer; even power adapters or accumulators of identical size and specification from other manufacturers are not permitted. Do not disassemble the laboratory device and its accessories (e.g., power adapters, cables, racks, accumulators or batteries) any further than described in the user manual!
- Only autoclave the laboratory device when this is permissible according to the user manual.

2.2 Purpose

This is an air displacement pipette for pipetting liquids of low to medium viscosity.

2.3 Limitations of use

The pipette is used for dispensing liquids within the following limits:

- Operating temperature from +15 °C to +40 °C (59 °F to 104 °F). Additional temperatures upon request
- vapor pressure up to 500 mbar

- Viscosity: 260 mPa s

With viscous liquids, the dispensing speed may need to be adjusted.

2.4 Application restrictions

Viscous and wetting liquids may compromise volumetric accuracy. Volumetric accuracy may also be affected when pipetting liquids whose temperature deviates from the ambient temperature by more than $\pm 1\text{ }^{\circ}\text{C}/\pm 1.8\text{ }^{\circ}\text{F}$.

2.5 Operating exclusions

The user must verify the suitability of the instrument for the intended purpose, as aggressive liquids and their vapors can damage the instrument (corrosion!). The instrument cannot be used for the following liquids:

- Liquids with very high steam pressure
- Liquids that corrode the following materials:
 - Fluoroelastomer rubber (FKM)
 - Polycarbonate (inspection window)
 - Polyetheretherketone (PEEK)
 - Polyoxyethylene (POM)
 - Polyphenylene sulfide (PPS) (with variable 50 μl instrument, 5 ml, 10 ml pipetting unit)
 - Polypropylene (PP)
 - Polyvinylidene fluoride (PVDF)

Additional information on the chemical resistance of plastics can be found at www.brand.de.

3 Functional and operational components

- 1** Pipetting button
3 Easy Calibration function
5 Pipetting shaft
7 Handle
9 Volume-change protection

- 2** Volume-setting wheel
4 Finger rest
6 Tip cone
8 Volume display
10 Tip ejection key

Label window

The instrument can be individually labeled on the finger rest:

- Remove the label window on the finger rest.
- Mark the labeling film.
- Reinsert the labeling film with window.


4 Pipetting

1. Inserting a tip

NOTICE

- > 2 ml, 5 ml and 10 ml devices should only be used with a built-in PE filter; see UV sterilization, p. 46.
- > Perfect analysis results can only be achieved by using quality tips. We recommend BRAND pipette tips. For additional information, refer to the accuracy table Accuracy table, p. 42.
- > Pipette tips are disposable products!

a.

a. Insert tips vertically:

Use the correct tips, in accordance with the volume range or color code!

Make sure that the tips are firmly in place and leak tight.

2. Setting the volume

a.


- a. Slide the volume-change protection upward (UNLOCK).
- b. Turn the volume-setting wheel to select the desired volume. In doing so, turn the adjustment wheel steadily, avoiding abrupt turning motions.
- c. Slide the volume-change protection downward (LOCK). The volume-setting wheel becomes noticeably more difficult to turn, but movement is not completely blocked.

3. Aspirating a sample

NOTICE

We recommend rinsing the pipette tip 5 times with the liquid before pipetting (aspirate and dispense liquid 5 times) to achieve maximum precision and accuracy. See also Checking the volume, p. 41.

- Press the pipetting button until first resistance is felt.
- Hold the device vertically and immerse the tip in the liquid.

- Allow the pipetting button to steadily move back to its original position.


Leave the tip immersed in the liquid for a few seconds, so that the set volume is aspirated completely. This is especially important when pipetting viscous media and when using pipettes with large volumes.

Volume range	Immersion depth [mm]	Wait time [s]
0.1 µl - 1 µl	1 - 2	1
1 µl - 100 µl	2 - 3	1
100 µl - 1000 µl	2 - 4	1
> 1,000 µl	3 - 6	3

NOTICE

Do not lay the instrument down when the tip is filled; this can cause the medium to flow into the instrument and contaminate it! The instrument should always be stored in the provided shelf mount or table stand and kept in an upright position, without any tip inserted.

4. Dispensing a sample

- Place the pipette tip against the vessel wall. Hold the pipette at an angle of 30-45° to the vessel wall.
- Press the pipetting button at a uniform speed until the first resistance is felt and hold it. To improve accuracy, comply with the corresponding wait time for serums, highly-viscous or low-density media.


- Empty the tip completely by overtravel: Press the pipetting button to the second stop.
- While doing this, wipe the pipette tip against the vessel wall.
- Remove the pipette tip from the vessel wall and allow the pipetting button to move back to its original position.

5. Ejecting a tip

NOTICE

The device should always be stored in the provided shelf mount or table stand and kept in an upright position, without any tip inserted.

- Hold the pipette shaft over a suitable disposal bin and press the tip ejection button all the way down.

5 Checking the volume

We recommend testing the instrument every 3 to 12 months depending on the level of use. However, the testing cycle can be adapted to meet individual requirements. The complete testing procedure (SOP) can be downloaded at www.brand.de.

You can download the detailed test instructions (SOP) at www.brand.de. For GLP- and ISO-compliant evaluations and documentation, we recommend the EASYCAL™ calibration software from BRAND. A demo version can be downloaded from <https://shop.brand.de/>.

Gravimetric volume testing of the pipette is carried out according to the following steps and complies with DIN EN ISO 8655:2022.

- a. Set the maximum specified instrument volume (for procedure, see Pipetting, p. 38).
- b. Condition the pipette before testing by aspirating and dispensing the test liquid (distilled water) with a pipette tip five times.
- c. Aspirate the test liquid and pipette into the weighing vessel.
- d. Weigh the pipetted amount with an analysis scale. (refer to the operating manual of the balance manufacturer.)
- e. Calculate the pipetted volume. In doing so, take into account the temperature of the test liquid.
- f. At least 10 pipetting series and weighings in three volume ranges (100%, 50%, 10%) are recommended. Two tips must be used for each volume range to be tested.

Calculation (for nominal volume)

x_i = weighing results

n = number of weighings

V_0 = nominal volume

Z = Correction factor (e.g. 1.0029 µl/mg at 20°C, 1013 hPa)

Mean:

$$\bar{x} = \frac{\sum x_i}{n}$$

Mean volume:

$$\bar{V} = \bar{x} * Z$$

Accuracy*:

$$A\% = \frac{\bar{V} - V_0}{V_0} * 100$$

Standard deviation*:

$$s = Z * \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

Coefficient of variation*:

$$CV\% = \frac{100 s}{\bar{V}}$$

*) Accuracy and coefficient of variation are calculated according to the formulas of statistical quality control.

6 Accuracy table

Transferette® S, adjustable

Volume range [µl]	Partial volume [µl]	A* ≤ ± %	CV* ≤ %	Sub steps [µl]	Recommended tip type [µl]
0.1 - 1	1	2	1.2	0.001	0.1 - 20
	0.5	4	2.4		
	0.1	20	12		
0.1 - 2.5	2.5	1.4	0.7	0.002	0.5 - 20
	1.25	2.5	1.5		
	0.25	12	6		
0.5-10	10	1	0.5	0.01	0.5 - 20
	5	1.6	1		
	1	7	4		
2 - 20	20	0.8	0.4	0.02	2 - 200
	10	1.2	0.7		
	2	5	2		
5 - 50	50	0.8	0.3	0.05	2 - 200
	25	1.2	0.5		
	5	4	2		
10 - 100	100	0.6	0.2	0.1	2 - 200
	50	0.8	0.4		
	10	3	1		
20 - 200	200	0.6	0.2	0.2	2 - 200
	100	0.8	0.3		
	20	3	0.6		
100 - 1000	1000	0.6	0.2	1	50 - 1000
	500	0.8	0.3		
	100	3	0.6		
500...5,000	5000	0.6	0.2	5	500...5,000
	2500	0.8	0.3		
	500	3	0.6		
1000 - 10,000	10000	0.6	0.2	10	1,000...10000
	5000	0.8	0.3		
	1000	3	0.6		

*A = Accuracy, CV = Coefficient of Variation

Transferpette® S, fixed-volume

Volume range [µl]	A* ≤ ± %	CV* ≤ %	Recommended tip type [µl]
10	1	0.5	0.5 - 20
20	0.8	0.4	2 - 200
25	0.8	0.4	2 - 200
50	0.8	0.4	2 - 200
100	0.6	0.2	2 - 200
200	0.6	0.2	2 - 200
500	0.6	0.2	50 - 1000
1000	0.6	0.2	50 - 1000
2000	0.8	0.3	500...5,000


*A = Accuracy, CV = Coefficient of Variation

Final test values based on the nominal volume (= max. volume) printed on the device and the specified partial volumes at the same temperature (20 °C/68 °F) of the device, surroundings and distilled water, in accordance with DIN EN ISO 8655.

7 Adjustment – Easy Calibration

The instrument is permanently calibrated for aqueous solutions. If it is determined that the pipette is operating inaccurately or to adjust the instrument to work with solutions of varying density and viscosity or with specially-shaped pipette tips, it can be calibrated using the Easy Calibration Technique.

b.

- Perform a volume check and determine the actual value; see Checking the volume.
- Remove label window (1) and labeling film (2): Gently move the clamp and lift it off.
- Using a paper clip or an unused pipette tip, remove the protective film (3) (the protective film can be discarded).

d.

- Slide the red adjustment slider back completely, lift the volume-setting wheel (decoupling) and release the adjustment slider.

e.

- Set the adjustment value:
Transferpette® S, adjustable: with the volume-setting wheel in the UNLOCK position, set to the previously determined actual value.
Transferpette® S, fixed-volume: set the volume by rotating in the +/- direction.
A volume check is recommended after every adjustment.

f.

- f. Slide the adjustment slider completely back again, push the volume-setting wheel downward and release the adjustment slider. Re-attach the labeling film and re-assemble the label window.

NOTICE

The change to factory settings is indicated by the red adjustment slider now visible in the label window.

8 Disinfection/autoclaving

8.1 Autoclaving

NOTICE Carry out a self-test of the effectiveness of the autoclave.

Maximum safety is achieved through vacuum sterilization. We recommend using sterilization bags.

- a. Eject the pipette tip.
- b. Fully unlock the volume adjustment lock (UNLOCK).
- c. Set the volume to the full value (e.g. 11.25 or 11.26) and leave the volume adjustment lock unlocked. This prevents jamming or damaging the volume setting during autoclaving.

Correct

Incorrect

- d. Pack the instrument in a sterilization bag, observing any applicable packaging regulations.
- e. Autoclave the complete pipette without further disassembly. Recommendation for autoclaving according to DIN EN 285 (see table below).
- f. Allow the pipette to completely cool and dry.

Temperature	121 °C
Pressure	2 bar
Holding time in autoclaves	15 min

After autoclaving, tighten the connection between the hand grip and the pipette shaft if necessary.

8.2 UV sterilization

The device is resistant to normal exposure to a UV disinfection lamp. The effects of the UV exposure may cause some color change.

8.3 PE filter

PE filter for Transferpette® S 2 ml, 5 ml + 10 ml:

A hydrophobic PE filter is used to prevent liquid from entering the pipette.

Change the filter if it becomes wet or contaminated.

- a. Use a flat object, such as a screwdriver.
- b. Remove the filter carefully, without damaging the tip cone.

Remove the filter before autoclaving!

The device can also be operated without a filter.

9 Maintenance

- a. Check the pipette tip cone for damage.
- b. Inspect the piston and seal for contamination.
- c. Check the device for leaks.


We recommend using the BRAND leak detector, the BRAND PLT unit. As an alternative to this, aspirate a sample and hold the device vertically for approx. 10 s. If a drop forms at the pipette tip, follow the malfunction remedy, see Troubleshooting, p. 50.

9.1 Check pipetting unit

- a. Check the pipette tip cone for damage.
- b. Inspect the piston and seal for contamination.
- c. Check the device for leaks.

We recommend using the BRAND leak detector, the BRAND PLT unit. As an alternative to this, aspirate a sample and hold the device vertically for approx. 10 s. If a drop forms at the pipette tip, follow the malfunction remedy, see Troubleshooting, p. 50.

9.2 Disassembling/cleaning (up to 1,000 µl)

- a. Detach the pipette shaft (S) from the hand grip by unscrewing it.
- b. Unscrew the upper part of the ejector unit (A) from the pipette shaft.
- c. Pull out the shaft (B, C and D) from the lower part of the ejector unit (E).
- d. Unscrew the piston unit (B).

NOTICE

Do not disassemble the piston unit (B) any further!

- e. Remove the spring (C) (not possible on Transferpette® S 1 µl, 2.5 µl and 10 µl!).
- f. Clean the parts shown with a soap solution or isopropanol, and then rinse with distilled water.
- g. Dry the parts (max. 120 °C/248 °F).
- h. Grease piston and seal with a very thin layer of supplied silicone grease.

Reassemble the cooled parts in reverse order. Only hand-tighten the piston unit and the upper part of the ejector unit (A, B).

9.3 Disassembly/cleaning (2–10 ml)

- Remove the entire shaft (S) from the hand grip by rotating at the upper end of the ejector (F) and remove the filter (K) from the bottom part of the shaft (H).
- Separate the bottom part of the ejector (F') by unscrewing it from the upper part of the ejector (F).
- Unscrew and dismantle the piston unit (G) with the ejector spring (I) and the bottom part of the shaft (H).
- Remove the O-ring-seal from the piston unit and clean it.

NOTICE

Do not disassemble the piston unit (G) any further!

- Clean the piston unit (G) and the bottom part of the shaft (H) with a soap solution or isopropanol, and then rinse with distilled water.
- Dry the parts (max. 120 °C/248 °F) and allow them to cool.
- Carefully lubricate the inside and outside of the O-ring (O) and mount it on the piston.

Reassemble the individual components in reverse order.

10 Troubleshooting

Problem	Possible cause	Corrective action
Tip dripping (device leaking)	Unsuitable tip	Only use high-quality tips
	Tip not seated tightly	Firmly press tip on
The instrument does not aspirate or aspirates too little; the dispensed volume is too low	Seal contaminated	Clean seal
	Seal or cone is damaged	Replace seal or shaft
	Piston is contaminated or damaged	Clean or replace piston
Aspiration is very slow	Shaft is clogged	Clean shaft
	Filter contaminated on 2 ml, 5 ml or 10 ml instruments	Change filter
Dispensed volume too large	Pipetting button pressed too far (to the over-stroke point) before aspirating	Ensure proper handling.
Piston sluggish	Piston is contaminated or not greased	Clean piston and apply grease

11 Product markings

Symbol or number	Meaning
	Read the user manual.
XXZXXXXX	Serial number
	The instrument is marked in accordance with the German Measurement and Calibration Act as well as the Measurement and Calibration Regulation. Character sequence DE-M (DE for Germany), framed by a rectangle, as well as the two last digits of the year the marking was affixed.
Data matrix code or quick response code	The codes link to the BRAND MyProduct website.
www.brand.de/ip	Hyperlink to BRAND patent site

12 Order Information

12.1 Order info/accessories

Transferpette® S, fixed-volume

Volume	Description	Order No.
10 µl	F-10	705808
20 µl	F-20	705816
25 µl	F-25	705820
50 µl	F-50	705828
100 µl	F-100	705838
200 µl	F-200	705844
500 µl	F-500	705854
1,000 µl	F-1000	705862
2000 µl	F-2000	705864

Transferpette® S, adjustable

Volume	Description	Order No.
0.1–1 µl	D-1	705868
0.1–2.5 µl	D-2.5	705869
0.5–10 µl	D-10	705870
2–20 µl	D-20	705872
5–50 µl	D-50	705873
10–100 µl	D-100	705874
20–200 µl	D-200	705878
100–1,000 µl	D-1000	705880
0.5–5 ml	D-5000	705882
1–10 ml	D-10000	705884


Accessories

	Description	Order No.
	Table stand for 6 single-channel or multi-channel pipettes.	704807
	Table stand for 1 single-channel or multi-channel pipette.	703440

	Description	Order No.
	Wall mount	704812
	Shelf/rack mount	704811

12.2 Spares

12.2.1 Transferpette® S up to 1,000 µl

Transferpette® S, fixed-volume


Volume	A	B	C	D	E
10 µl	705508	704601	-	704721*	704739
20 µl	705509	704602	704610	704723	704740
25 µl	705509	704608	704614	704723	704741
50 µl	705509	704654	704661	704724	704742
100 µl	705509	704654	704661	704724	704743
200 µl	705509	704655	704662	704725	704745
500 µl	705511	704656	704663	704726	704746
1000 µl	705511	704656	704663	704726	704747

Transferpette® S, adjustable

Volume	A	B	C	D	E
0.1– 1 µl	705508	704600	-	704718*	704730
0.1– 2.5 µl	705508	704616	-	704719*	704731
0.5– 10 µl	705508	704601	-	704721*	704732
2– 20 µl	705509	704602	704610	704723	704733
5– 50 µl	705509	704615	704617	704722	704734
10– 100 µl	705509	704654	704661	704724	704735
20– 200 µl	705509	704655	704662	704725	704736
100– 1000 µl	705511	704656	704663	704726	704737

* Seal permanently installed in shaft – not removable!

12.2.2 Transferpette® S, 2 ml, 5 ml and 10 ml

NOTICE

The appearance and dimensions of the spare parts correspond to the respective nominal volume.

Transferpette® S, fixed-volume

Volume	F + F'	G	H	I	O
2 ml	704765	704606	703247	704626	7288

Transferpette® S, adjustable

Volume	F + F'	G	H	I	O
0.5– 5 ml	704766	704606	703247	704626	7288
1– 10 ml	704767	704607	704628	704626	7298

12.3 Additional accessories

Description	Order No.
Filter for 2 ml and 5 ml volumes, PU 25 pc.	704652
Filter for 10 ml volume, PU 25 pc.	704653
Silicone grease for volumes up to 1000 µl	705502
Silicone grease for volumes 2 ml / 5 ml / 10 ml	703677
Label window, PU 1 pc.	704750
Labeling film, PU 5 pc.	704751
PLT unit (pipette leak detector)	703970

13 Repairs

13.1 Sending for repair

NOTICE

Transporting hazardous materials without approval is prohibited by law.

Clean the instrument thoroughly and decontaminate!

- When returning products, please enclose a general description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

Outside USA and Canada

Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de.

Within USA and Canada

Please clarify the requirements for the return delivery with BrandTech Scientific, Inc **before** sending the instrument in for service.

Return only cleaned and decontaminated instruments to the address provided with the Return Authorization Number. Place the Return Authorization number so that it is clearly visible on the outside of the package.

Contact addresses

Germany:

BRAND GMBH + CO KG
Otto-Schott-Straße 25
97877 Wertheim (Germany)
T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

USA and Canada:

BrandTech® Scientific, Inc.
11 Bokum Road
Essex, CT 06426-1506 (USA)
T +1-860-767 2562
F +1 - 860 - 767 2563
info@brandtech.com
www.brandtech.com

India:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park,
Powai
Mumbai-400 076 (India)
T +91 22 42957790
F +91 22 42957791
info@brand.co.in
www.brand.co.in

China:

BRAND (Shanghai) Trading Co., Ltd.
Rm 201-202, North Tower,
No. 199 Kaibin Rd, Xuhui District, Shanghai
Shanghai 200030 (P.R. China)
T +86 21 6422 2318
F +86 21 6422 2268
info@brand.com.cn
www.brand.cn.com

14 Calibration service

The ISO 9001 and GLP guidelines require regular inspection of your volume measuring devices. We recommend performing a volume check every 3 to 12 months. The cycle is dependent on the individual requirements of the device. Checks should be performed more frequently, in case of high frequency of use or the use of aggressive Liquids.

The complete SOP for testing can be downloaded from www.brand.de or www.brandtech.com.

BRAND also offers you the option of having your devices calibrated through our factory calibration service or through our accredited calibration laboratory. Just send us the devices to be calibrated, indicating the type of calibration you would like. You will get your devices back in a few days. A detailed calibration report (factory calibration) or an accredited calibration certificate in accordance with DIN EN ISO/IEC 17025 is enclosed with each device. More information can be obtained from your retailer or directly from BRAND. The order document is available for download at www.brand.de (Service & Support).

For customers outside Germany

If you would like to use our calibration service, please contact one of our service partners in your region. Our service partners can forward your devices to BRAND for factory calibration, if required.

15 Information about your laboratory instrument

The online service MyProduct (<https://www.brand.de/myproduct>) offers quality certificates, accessories, and technical documentation for your laboratory instrument Transferpette® S. By entering the serial and article number there, you will receive information specific to your individual instrument.

On Transferpette® S, you will also find serialized information encoded in the quick response code. This also links to the MyProduct website, where you can access this operating manual and additional certificates for your instrument.

You will also find a data matrix code on some instruments (Transferpette® S, HandyStep® touch, and HandyStep touch® S). Scan this using a standard reader app to access the above-mentioned information via the URL <https://www.brand.de/myproduct>.

16 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the instrument or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

USA and Canada:

Find more warranty information on www.brandtech.com.

17 Disposal

Before disposal, observe the relevant national disposal regulations and ensure that the product is disposed of properly.

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