[Kurzanleitung](#) | [Short instructions](#) | [Instruction brève](#) |
[Guía rápida](#) | [Guida rapida](#) | [Breves instruções](#) | [快速指南](#)

Transferpette® pro

Mikroliterpipetten | Micropipettes

Impressum

BRAND GMBH + CO KG

Otto-Schott-Str. 25
97877 Wertheim (Germany)

T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

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Table of contents

1	Use operating manual	33
1.1	Signal words and their meaning	33
1.2	Presentation of descriptions of actions	33
2	Safety regulations.....	34
2.1	Safety regulations	34
2.2	Target group.....	35
2.3	Purpose	36
2.4	Use	36
2.5	Limitations of use.....	36
2.6	Operating exclusions	36
2.7	Improper use	40
2.8	Foreseeable misuse	40
3	Functional and operational components.....	41
4	Functional and operational components.....	45
5	Pipetting	49
6	Pipetting	53
7	Adjustment.....	57
8	Order Information.....	59

1 Use operating manual

- Please carefully read the operating manual before using the device for the first time.
- Keep the operating manual in an easily accessible place. It is part of the instrument.
- Be sure to include the operating manual if you transfer possession of this device to a third party.
- You can find updated versions of the instructions for use at www.brand.de/om.
- You can find updated versions of the Instructions for Use at <https://www.vitlab.com/manuals>.

1.1 Signal words and their meaning

Signal word	Meaning
⚠ WARNING or ⚠WARNING! ...	WARNING indicates a dangerous situation that, if not avoided, could result in death or serious injury.
⚠ CAUTION or ⚠CAUTION! ...	CAUTION indicates a hazardous situation that, if not avoided, could result in moderate or minor injury.
NOTICE or NOTICE ...	NOTE is used to address actions that are not related to physical injury. Example: Possible property damage.

1.2 Presentation of descriptions of actions

Format	Meaning
1. Task	Indicates a task.
a., b., c.	Indicates the individual steps of a task.
>	Indicates a prerequisite for a task.
⇒	Indicates a result of a completed task.

2 Safety regulations

2.1 Safety regulations

Please read carefully!

The laboratory device Transferpette® pro Can be used in combination with hazardous materials, operations and equipment. However, the instructions for use cannot identify all the safety problems that may occur. It is the responsibility of the user to ensure compliance with safety and health regulations and to determine the appropriate restrictions before use.

- Each user must have read and observe the instructions for use accompanying the laboratory device before using the device. The laboratory equipment must only be used by trained and qualified personnel.
- Follow general hazard warnings and safety regulations, e.g. wear protective clothing, eye protection and protective gloves.
- When working with infectious or hazardous samples/media (e.g. hazardous materials), the general safety rules in the laboratory must be followed and regulations for handling the samples/media must be observed. The information provided by the media manufacturers (e.g. safety data sheets) must be observed.
- Use the laboratory equipment only for pipetting or dispensing media within the defined limits and limitations of use. Observe usage exclusions.
- When working with flammable media, take precautions to avoid electrostatic charges, e.g. do not dose into plastic containers and do not rub equipment with a dry cloth. Do not use the laboratory equipment in explosive atmospheres. If in doubt, it is essential to contact the manufacturer or dealer.

- Always check that the laboratory equipment is in good condition before use. If there are any signs of malfunctions in the laboratory equipment (e.g. Stiff piston, leaks or power supply), immediately stop working with the appliance and refer to the troubleshooting section in the user manual. If necessary, contact the manufacturer.
- Always work in such a way that neither users nor other persons are at risk. Avoid splashes. Use only suitable vessels. Never use unnecessary force or force while operating, cleaning, or maintaining the laboratory equipment.
- If the laboratory equipment is powered by a power supply, batteries or rechargeable batteries, the correct condition of the components and the connection to the equipment must be checked regularly. Do not operate the laboratory equipment or its accessories in an unprotected, damp or wet environment.
- Do not make any technical changes. Use only original manufacturer's spare parts, and do not use power supplies or batteries of identical sizes and specifications from other manufacturers. Do not disassemble the laboratory equipment and its accessories (e.g., power supplies, cables, stands, batteries, etc.) any further than the instructions for use.
- Do not autoclave the laboratory equipment unless permitted by the instructions for use.

2.2 Target group

The operating manual is intended for users who use the laboratory instrument in the course of their professional activities. Users are familiar with the typical safety regulations and working methods in laboratories and have been trained accordingly. They can recognize potential hazards and protect themselves from them. The operating manual assumes this expertise and does not replace basic laboratory training or specific safety training.

2.3 Purpose

This is an air displacement pipette for pipetting liquids of medium density and low to medium viscosity.

2.4 Use

Use the laboratory device Transferpette® pro only for pipetting or dispensing liquids within the defined limits of use.

2.5 Limitations of use

The pipette is used for dispensing liquids within the following limits:

- Operating temperature from +15 °C to +40 °C (59 °F to 104 °F).
Additional temperatures upon request
- vapor pressure up to 500 mbar
- Viscosity: 260 mPa s

With viscous liquids, the dispensing speed may need to be adjusted.

2.6 Operating exclusions

The user must verify the suitability of the instrument for the intended purpose, as aggressive liquids and their vapors can damage the instrument (corrosion!). The instrument cannot be used for the following liquids:

- Liquids with very high steam pressure
- Liquids that corrode the following materials:
 - Fluoroelastomer rubber (FKM)
 - Polycarbonate (inspection window)
 - Polyetheretherketone (PEEK)
 - Polyoxymethylene (POM)
 - Polyphenylene sulfide (PPS) (with variable 50 µl instrument, 5 ml, 10 ml pipetting unit)
 - Polypropylene (PP)
 - Polyvinylidene fluoride (PVDF)

Additional information on the chemical resistance of plastics can be found at www.brand.de.

- Liquids with very high steam pressure
- Liquids that corrode the following materials
 - Fluoroelastomer rubber (FKM)
 - Polycarbonate (PC)
 - Polyphenylsulphide (PVDF) (with 50 µl, 100 µl, 200 µl and 300 µl instruments)
 - Polypropylene (PP)
 - Polyvinylidene fluoride (PVDF)
 - Silicone

Additional information on the chemical resistance of plastics can be found at www.brand.de.

- Very high vapor pressure fluids
- Liquids that attack the following materials:
 - Fluoro-elastomer rubber (FKM)
 - Polycarbonate (PC) (viewing window)
 - Polyetheretherketone (PEEK)
 - Polyphenylene sulfide (PPS) (50 µl pipette)
 - Polypropylene (PP)
- Liquids with very high steam pressure
- Liquids that corrode the following materials:
 - Polypropylene (PP)
 - Polycarbonate (PC) (inspection window)
 - Fluoroelastomer rubber (FKM)
 - Polyetheretherketone (PEEK)
- Liquids with very high steam pressure
- Liquids that corrode the following materials:
 - Polycarbonate (PC)
 - Polyphenylene sulfide (PPS)
 - Polyvinylidene fluoride (PVDF)
 - Silicone

- Liquid with very high steam pressure
- Liquids that corrode the following materials
 - Fluoroelastomer rubber (FKM)
 - Polycarbonate/polybutylene terephthalate (PC/PBT)
 - Polyetheretherketone (PEEK)
 - Polypropylene (PP)
 - Polyvinylidene fluoride (PVDF)

NOTICE

Autoclaving

The hand grip is not autoclavable!

- > Heat and steam during autoclaving will damage the electronics.
-

Additional information on the chemical resistance of plastics can be found at www.brand.de.

- Liquids with very high steam pressure
- Liquids that corrode the following materials:
 - Fluoroelastomer rubber (FKM),
 - Polycarbonate (PC)
 - Polycarbonate/polybutylene terephthalate (PC/PBT)
 - Polyetherimide (PEI)
 - Polypropylene (PP)
 - Polyvinylidene fluoride (PVDF)
 - Silicone

NOTICE

Autoclaving

The hand grip is not autoclavable!

- > Heat and steam during autoclaving will damage the electronics.
-

Additional information on the chemical resistance of plastics can be found at www.brand.de.

- Liquids with very high steam pressure
- Liquids that corrode the following materials:
 - Fluoroelastomer rubber (FKM)

Polycarbonate (PC) (inspection window)

Polyetheretherketone (PEEK)

Polypropylene (PP)

Polyvinylidene fluoride (PVDF)

- for liquids with very high steam pressure
- Liquids that corrode the following materials:

Fluoroelastomer rubber (FKM)

Polyamide (PA)

Polycarbonate (inspection window)

Polyetheretherketone (PEEK)

Polyphenylene sulfide (PPS)

Polypropylene (PP)

Polyvinylidene fluoride (PVDF)

Additional information on the chemical resistance of plastics can be found at www.brand.de.

- for liquids with very high steam pressure
- Liquids that corrode the following materials:

Fluoroelastomer rubber (FKM)

Polyamide (PA)

Polycarbonate (inspection window)

Polyetheretherketone (PEEK)

Polyphenylene sulfide (PPS)

Polypropylene (PP)

Polyvinylidene fluoride (PVDF)

Silicone (SI)

Additional information on the chemical resistance of plastics can be found at www.brand.de.

For more information on chemical resistance of plastics, visit www.brand.de.

2.7 Improper use

Various risks may arise if the laboratory instrument is used improperly. These risks include: inaccurate liquid delivery, damage to the laboratory instrument, and the risk of contamination, infection, and injury from contact with the pipetted media.


Any use other than for pipetting or dispensing liquids within the defined operating limits is considered improper use.

2.8 Foreseeable misuse

A typical misuse is pipetting or dispensing liquids with too high viscosity or using unsuitable tips.

3 Functional and operational components

Front

Adjustable finger rest	The laboratory instrument Transferpette® pro has an adjustable finger rest. This allows you to adapt the pipette to your preferred grip, see "Pipetting The instrument can be labeled at the finger rest: to do so, remove the label window from the finger rest and take out the label strip.
Volume-change protection	The volume adjustment lock prevents the volume from being changed during pipetting work. To unlock it, slide the volume adjustment lock toward the pipetting button.
Volume display	The numbers in the display are read from top to bottom; the white dash corresponds to the decimal point.
Counter mechanism protection	Once the volume adjustment lock is released, set the volume using the volume-setting wheel. The integrated counter mechanism protection overrides the volume-setting wheel once the maximum or minimum volume setting is reached: the volume-setting wheel remains rotatable but no longer adjusts the volume.
Hand grip	Screw the pipetting unit into the hand grip. Attach the tip onto the tip cone.

Rear

Permanent adjustment to factory settings:
Easy Calibration

Temporary adjustment to
changing conditions:
User Adjustment

User Adjustment scale

Cover

Seal

Label window

Serial number and
product markings

QR code:

Use your smartphone, tablet,
or webcam to scan and visit
the following website:
www.brand.de/myproduct

The website contains
serialized information about
your laboratory instrument.

Easy Calibration Technology

The Easy Calibration function is located in the finger rest and is used to reset the pipette to its factory settings (see Adjustment, p. 57).

User Adjustment Technology

The hand grip also includes the User Adjustment function. This allows the pipette to be adjusted for specific liquids and dispensing conditions. See Temporary adjustment: User Adjustment

The switch for user adjustment is located behind the cover. A seal is applied at delivery. Remove it upon first use and dispose of it.


QR Code and serialized information

The QR code provides access to www.brand.de/myproduct and serialized information about your pipette.

If you wish to access the MyProduct information without the QR code, you will also need the order number and serial number of your pipette.

4 Functional and operational components

Front

Adjustable finger rest

The laboratory instrument Transferpette® pro has an adjustable finger rest. This allows you to adapt the pipette to your preferred grip, see "Pipetting

The instrument can be labeled at the finger rest: to do so, remove the label window from the finger rest and take out the label strip.

Volume-change protection

The volume adjustment lock prevents the volume from being changed during pipetting work. To unlock it, slide the volume adjustment lock toward the pipetting button.

Volume display

The numbers in the display are read from top to bottom; the white dash corresponds to the decimal point.

Counter mechanism protection

Once the volume adjustment lock is released, set the volume using the volume-setting wheel. The integrated counter mechanism protection overrides the volume-setting wheel once the maximum or minimum volume setting is reached: the volume-setting wheel remains rotatable but no longer adjusts the volume.

Hand grip

Screw the pipetting unit into the hand grip. Attach the tip onto the tip cone.

Rear

Permanent adjustment to factory settings:
Easy Calibration

Temporary adjustment to
changing conditions:
User Adjustment

Label window

Step value scale

Serial number and
product markings

Cover

Seal

QR code:

Use your smartphone, tablet,
or webcam to scan and visit
the following website:
www.brand.de/myproduct

The website contains
serialized information about
your laboratory instrument.

Easy Calibration Technology

The Easy Calibration function is located in the finger rest and is used to reset the pipette to its factory settings (see Adjustment, p. 57).

User Adjustment Technology

The hand grip also includes the User Adjustment function. This allows the pipette to be adjusted for specific liquids and dispensing conditions. See Temporary adjustment: User Adjustment

The switch for user adjustment is located behind the cover. A seal is applied at delivery. Remove it upon first use and dispose of it.

QR Code and serialized information

The QR code provides access to www.brand.de/myproduct and serialized information about your pipette.

If you wish to access the MyProduct information without the QR code, you will also need the order number and serial number of your pipette.

5 Pipetting

a.

- a. Turn the finger rest to a comfortable working position.

b.

- b. Slide the volume adjustment lock in the direction shown, against slight resistance.

c.

- c. Set the volume using the volume-setting wheel.

NOTICE If the volume-setting wheel is turned beyond the maximum or minimum volume, it glides over the volume adjustment and thus protects the counter mechanism from damage.

d.

d. Close the volume adjustment lock.

⇒ The volume-setting wheel can still be turned but will no longer change the volume.

e.

e. Attach the pipette tip. Ensure a secure fit.

The 2 ml to 10 ml pipettes should only be used with a built-in PE filter (see UV sterilization).

NOTICE Pipette tips are disposable products!


f.

f. Press the pipetting button down to the first stop.

NOTICE We recommend rinsing the pipette tip five times with the liquid (aspirate and dispense) before pipetting in order to ensure maximum precision and accuracy.

g.

g. Immerse the tip in the liquid.

h.


h. Slowly release the pipetting button.

⇒ Liquid is aspirated.

Keep the tip submerged until the volume is fully aspirated. Increase the wait time when pipetting more viscous liquids or larger volumes.

⚠CAUTION! Do not lay the pipette down with filled tips. Contamination may occur!

Volume range	Immersion depth [mm]	Wait time [s]
0.1–1 µl	1–2	1
1–100 µl	2–3	1
100–1000 µl	2–4	1
> 1,000 µl	3–6	3

- i. To dispense the liquid, hold the tip at an angle against the vessel wall, slowly press the pipetting button, and wipe the tip.

To improve accuracy, comply with the corresponding wait time for serums, highly-viscous or low-density fluids.

To fully empty the tip, press the pipetting button down to the second stop (Fig. f).

- j. To remove the tip, hold the pipette over a container and press the tip ejection button.

Storage

You can also hook the Transferpette® pro in the holder or rack with an adjustable finger rest.

⚠CAUTION! Do not hang the pipette with a filled tip in the holder. Contamination may occur!

6 Pipetting

a.

- a. Turn the finger rest to a comfortable working position.

b.

- b. Slide the volume adjustment lock in the direction shown, against slight resistance.

c.

- c. Set the volume using the volume-setting wheel.

NOTICE If the volume-setting wheel is turned beyond the maximum or minimum volume, the counter mechanism protection prevents damage to the instrument: you can turn the volume-setting wheel, but it no longer changes the volume.

d.

- d. Close the volume adjustment lock.
NOTICE You can turn the volume-setting wheel, but it no longer adjusts the volume.

e.

- e. Attach the pipette tips. Ensure a secure fit.

The 2 ml to 10 ml pipettes should only be used with a built-in PE filter (see UV sterilization).


NOTICE Pipette tips are disposable products!

f.

- f. Press the pipetting button down to the first stop.

g.

g. Immerse the tips in the liquid.

h.

h. Slowly release the pipetting button.


⇒ Liquid is aspirated.

NOTICE The ISO 8655 standard requires that pipette tips are pre-wetted once before the actual pipetting procedure.

⚠CAUTION! Do not lay the pipette down with filled tips. Contamination may occur!

With more viscous liquids and larger volumes, keep the tip immersed for a waiting time to ensure the volume is fully aspirated.

Volume range	Immersion depth [mm]	Wait time [s]
0.1–1 µl	1–2	1
1–100 µl	2–3	1
100–1,000 µl	2–4	1
> 1,000 µl	3–6	3

- i. To dispense the liquid, hold the tips at an angle against the vessel wall, slowly press the pipetting button, and wipe the tips.

To improve accuracy, comply with the corresponding wait time for serums, highly-viscous or low-density fluids.

To fully empty the tips, press the pipetting button down to the second stop (Fig. f).

- j. To remove the tips, hold the pipette over a container and press the tip ejection button.

Storage


You can also hook the Transferpette® pro in the holder or rack with an adjustable finger rest.

⚠CAUTION! Do not hang the pipette with a filled tip in the holder. Contamination may occur!

7 Adjustment

The device is permanently calibrated for aqueous solutions. The pipette can be calibrated if it is operating inaccurately or must be adjusted to work with solutions of varying density and viscosity or with specially-shaped pipette tips.

b.

- a. Perform a volume check and determine the actual value; see Checking the volume.
- b. Removing the cover: Push the clamp forward, gently lift and pull backward.
- c. Using a paper clip or an unused pipette tip, remove the protective film (the protective film can be discarded).
- d. Slide the red adjustment slider back completely, lift the volume-setting wheel (decoupling) and release the adjustment slider.

d.

e.

- e. Set the adjustment value:
Make corrections by turning the volume-setting wheel in +/- direction. A volume check is recommended after every adjustment.

q.

- f. Slide the adjustment slider completely back again, push the volume-setting wheel downward and release the adjustment slider. Remount the cover

NOTICE

The change to factory settings is indicated by the red adjustment slider.

8 Order Information

Volume	Description	Cat. No.
0.1 - 2.5 µl	D-2.5	704969
0.5 - 10 µl	D-10	704970
2 - 20 µl	D-20	704972
5 - 50 µl	D-50	704973
10 - 100 µl	D-100	704974
20 - 200 µl	D-200	704978
100 - 1000 µl	D-1000	704980
0.5-5 ml	D-5000	704982
1-10 ml	D-10000	704984

Spares up to 1000 µl

For illustrations of spare parts B ... D, see Disassembling/cleaning (Transferpette® α up to 1000 µl)

Volume	B	C	D
0.1 - 2.5 µl	704616	–	704618*
0.5-10 µl	704601	–	704621*
2-20 µl	704602	704610	704622
5 - 50 µl	704615	704617	704659
10 - 100 µl	704654	704661	704623
20 - 200 µl	704655	704662	704624
100 - 1000 µl	704656	704663	704625

* 0.1 - 2.5 µl / 0.5-10 µl incl. seal

Replacement parts 5 ml/10 ml

For illustrations of spares G ... K, see Disassembling/cleaning (Transferpette® α 5 ml / 10 ml) .

Volume	G	H	I	K
0.5-5 ml	704606	703247	704626	704652*

Volume	G	H	I	K
1–10 ml	704607	704628	704626	704653*

* PU = 25 pcs.

Accessories

Description	Cat. No.
Shelf mount for 1 Transferpette® α	704809
Table stand for 6 Transferpette® α's	704806
Silicone grease for Transferpette® α up to 1,000 µl	705502
Silicone grease for Transferpette® α 5 ml/10 ml	703677

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