

 PRODUCT
INFORMATION

Liquid Handling Station

BRAND Liquid Handling Station


Contamination-free, automatic pipetting into 384-well microplates with the BRAND Liquid Handling Station

Introduction

With the BRAND Liquid Handling Station (LHS), a wide variety of everyday laboratory tasks can be performed automatically, with accurate results. Manually pipetting into 384-well microplates, as is common in high-throughput screening (HTS), is easily susceptible to errors and contamination. The objective of this technical note is to determine the contamination risk for pipetting into a 384-well microplate with the LHS. A common method used to establish this risk is to pipette a measuring solution into a microplate in a checkerboard pattern. Because the measuring solution is missing in every second well, contamination between wells can be measured.

Material and method

Every second well of a 384-well microplate (Cat. No. 781620), beginning with Well A2, is filled in checkerboard pattern with 50 μ l of a 1.5% solution of the food coloring Patent Blue V (E131). Next, each empty well is filled with 50 μ l of fully deionized (DI) water. The pipetting procedure is carried out once with the 50 μ l single-channel liquid end (SC LE) and once with the 50 μ l multi-channel liquid end (MC LE). In addition, in another plate, each well is filled with 50 μ l of DI water. This plate serves as the negative control. The absorption of all plates is then measured in the Nanoquant Infinite M200 Pro from Tecan at 640 nm. Based on the signal strengths in the wells, in which only DI water had been pipetted, it can then be determined if contamination of the food coloring has taken place.



384-well microplate (781620)

Results and discussion

Analysis of the negative control shows a mean absorption value of $\bar{\varnothing}_{Neg} = 0.0381 \pm 0.0013$. Because the measured values of the negative control are not normally distributed, the critical value, i.e. the value at which contamination can be assumed, is defined as the maximum measured value of the negative control. This value amounts to $A_{crit} = 0.0468$. This absorption corresponds to a substance quantity of 10.24 pmol and a coloring solution volume of $2.51 \cdot 10^{-3}$ pl.

As shown in Figure 1, no contamination was found in any of the wells filled with water in both 384-well microplates. For the pipetting procedure with the 50 µl SC LE, a maximum absorption of $A_{Max SC} = 0.0428$ was measured, and for the 50 µl MC LE, a maximum absorption of $A_{Max MC} = 0.0409$ was measured.

↔	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A	0,0373	3,1538	0,0378	3,1695	0,0370	3,0986	0,0400	3,1431	0,0375	3,1559	0,0376	3,0996	0,0375	3,1405	0,0377	3,1252	0,0368	3,1132	0,0373	3,1454	0,0377	3,0963	0,0378	3,0087
B	3,1366	0,0378	3,1083	0,0368	3,0732	0,0377	3,1206	0,0377	3,0990	0,0382	3,0976	0,0376	3,0569	0,0367	3,1512	0,0375	3,1327	0,0398	3,1388	0,0365	3,0804	0,0369	3,1498	0,0377
C	0,0374	3,0646	0,0381	3,1636	0,0375	3,0960	0,0378	3,1364	0,0380	3,0787	0,0379	3,1527	0,0403	3,0517	0,0377	3,1334	0,0369	3,1208	0,0372	3,1327	0,0369	3,0960	0,0375	3,1512
D	3,0520	0,0385	3,1204	0,0378	3,0926	0,0386	3,1329	0,0380	3,1064	0,0383	3,1634	0,0381	3,1214	0,0371	3,0688	0,0375	3,0763	0,0370	3,1355	0,0377	3,0902	0,0376	3,1281	0,0384
E	0,0382	3,0704	0,0380	3,0804	0,0376	3,1346	0,0388	3,1626	0,0379	3,1182	0,0388	3,1188	0,0373	3,0782	0,0376	3,0899	0,0373	3,0755	0,0370	3,1030	0,0374	3,0760	0,0379	3,1162
F	3,0636	0,0387	3,1943	0,0382	3,1162	0,0393	3,1752	0,0389	3,1019	0,0381	3,1157	0,0381	3,0797	0,0378	3,1463	0,0379	3,1477	0,0374	3,1431	0,0376	3,1149	0,0380	3,1277	0,0381
G	0,0396	3,0550	0,0392	3,1020	0,0393	3,1149	0,0396	3,1145	0,0390	3,1165	0,0386	3,1244	0,0385	3,0497	0,0377	3,1448	0,0379	3,0899	0,0377	3,1118	0,0381	3,0878	0,0383	3,1431
H	3,0716	0,0387	3,1284	0,0384	3,1286	0,0393	3,1117	0,0395	3,1112	0,0396	3,1581	0,0397	3,1161	0,0396	3,1156	0,0394	3,0614	0,0395	3,1265	0,0389	3,1700	0,0390	3,1202	0,0394
I	0,0411	3,0690	0,0394	3,1392	0,0401	3,1291	0,0399	3,1248	0,0396	3,0846	0,0394	3,1127	0,0394	3,0938	0,0377	3,2047	0,0391	3,1229	0,0403	3,1277	0,0388	3,0750	0,0394	3,1534
J	3,0484	0,0395	3,1570	0,0387	3,1326	0,0389	3,1437	0,0390	3,0228	0,0382	3,1419	0,0389	3,1111	0,0381	3,1875	0,0388	3,1243	0,0390	3,1539	0,0390	3,1413	0,0395	3,1223	0,0396
K	0,0386	3,0773	0,0390	3,0981	0,0379	3,1536	0,0381	3,1393	0,0406	3,1266	0,0383	3,1420	0,0378	3,0279	0,0381	3,1426	0,0387	3,0687	0,0395	3,1219	0,0384	3,0971	0,0390	3,1127
L	3,0830	0,0390	3,1248	0,0387	3,0816	0,0382	3,1569	0,0382	3,1255	0,0384	3,0374	0,0382	3,0695	0,0389	3,1395	0,0384	3,0905	0,0392	3,1219	0,0397	3,0360	0,0382	3,1044	0,0403
M	0,0388	3,0948	0,0386	3,1854	0,0397	3,0760	0,0380	3,1541	0,0389	3,1123	0,0381	3,1406	0,0391	3,0942	0,0402	3,1702	0,0388	3,1107	0,0394	3,1226	0,0393	3,0930	0,0383	3,1251
N	3,0809	0,0394	3,1259	0,0387	3,1104	0,0415	3,1007	0,0386	3,0967	0,0388	3,1230	0,0384	3,1248	0,0392	3,1101	0,0400	3,1086	0,0385	3,0812	0,0403	3,1493	0,0399	3,1197	0,0387
O	0,0395	3,0899	0,0389	3,1535	0,0380	3,0699	0,0389	3,1344	0,0392	3,1365	0,0385	3,0711	0,0393	3,0875	0,0393	3,1297	0,0385	3,0748	0,0428	3,1092	0,0390	3,0956	0,0411	3,1517
P	3,0064	0,0396	3,0461	0,0394	3,0673	0,0389	3,0818	0,0390	3,0744	0,0395	3,0662	0,0392	3,0455	0,0398	3,1341	0,0409	3,0838	0,0400	3,0942	0,0397	3,0782	0,0401	3,1485	0,0390

↔	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A	0,0372	3,2006	0,0376	3,2188	0,0367	3,2030	0,0403	3,1748	0,0376	3,1283	0,0372	3,1563	0,0371	3,1735	0,0371	3,1699	0,0370	3,2204	0,0377	3,2213	0,0373	3,1558	0,0381	3,1401
B	3,1520	0,0377	3,2011	0,0369	3,1442	0,0375	3,2353	0,0373	3,1557	0,0373	3,1574	0,0373	3,1394	0,0370	3,2194	0,0376	3,2011	0,0400	3,1711	0,0368	3,1528	0,0377	3,1668	0,0409
C	0,0369	3,1219	0,0373	3,1818	0,0366	3,1544	0,0369	3,2336	0,0375	3,2001	0,0374	3,0956	0,0360	3,1689	0,0367	3,2012	0,0366	3,1517	0,0381	3,1542	0,0369	3,1161	0,0378	3,1992
D	3,1398	0,0374	3,1806	0,0372	3,1963	0,0387	3,2714	0,0372	3,2526	0,0379	3,1805	0,0377	3,1081	0,0368	3,1578	0,0372	3,1365	0,0367	3,2151	0,0375	3,1216	0,0375	3,1590	0,0382
E	0,0368	3,1300	0,0369	3,1796	0,0361	3,2132	0,0374	3,2504	0,0371	3,1814	0,0373	3,1653	0,0363	3,1482	0,0369	3,1697	0,0361	3,1716	0,0361	3,2111	0,0367	3,1977	0,0378	3,1958
F	3,1490	0,0373	3,2118	0,0371	3,2291	0,0376	3,2496	0,0375	3,2466	0,0370	3,2299	0,0375	3,1849	0,0368	3,2113	0,0368	3,1847	0,0366	3,1973	0,0367	3,2130	0,0380	3,1117	0,0381
G	0,0376	3,1142	0,0372	3,2287	0,0378	3,1612	0,0379	3,2095	0,0375	3,1601	0,0372	3,1301	0,0367	3,1464	0,0367	3,1926	0,0366	3,1617	0,0367	3,1755	0,0370	3,1590	0,0389	3,1933
H	3,1456	0,0368	3,1920	0,0370	3,1749	0,0376	3,2268	0,0375	3,1583	0,0378	3,1946	0,0383	3,1313	0,0377	3,2067	0,0376	3,1627	0,0376	3,1653	0,0375	3,1752	0,0382	3,1256	0,0396
I	0,0384	3,2011	0,0388	3,2617	0,0388	3,1607	0,0376	3,1785	0,0379	3,1771	0,0378	3,1925	0,0374	3,1598	0,0364	3,1457	0,0380	3,1752	0,0378	3,1761	0,0374	3,1342	0,0385	3,2082
J	3,1566	0,0378	3,2251	0,0368	3,2255	0,0370	3,1915	0,0373	3,1571	0,0374	3,2442	0,0376	3,1418	0,0374	3,1476	0,0381	3,2094	0,0373	3,2310	0,0372	3,1561	0,0366	3,2032	0,0388
K	0,0371	3,1283	0,0371	3,2440	0,0362	3,1891	0,0368	3,1442	0,0378	3,1875	0,0376	3,1903	0,0361	3,1001	0,0369	3,1742	0,0374	3,2072	0,0396	3,2247	0,0368	3,1304	0,0372	3,2253
L	3,1414	0,0371	3,2274	0,0371	3,1866	0,0366	3,2300	0,0367	3,2093	0,0371	3,1624	0,0365	3,1603	0,0378	3,2407	0,0370	3,1895	0,0377	3,1733	0,0389	3,1876	0,0373	3,1834	0,0388
M	0,0368	3,1833	0,0364	3,1720	0,0375	3,1726	0,0366	3,2068	0,0370	3,2093	0,0366	3,1273	0,0373	3,0997	0,0381	3,2456	0,0368	3,1906	0,0375	3,1921	0,0371	3,2067	0,0367	3,2008
N	3,2058	0,0369	3,1870	0,0368	3,1882	0,0394	3,2202	0,0376	3,1561	0,0370	3,2052	0,0368	3,1405	0,0371	3,2287	0,0378	3,1419	0,0374	3,2603	0,0379	3,1701	0,0383	3,2031	0,0374
O	0,0370	3,2001	0,0369	3,2590	0,0357	3,1730	0,0361	3,1742	0,0372	3,2121	0,0365	3,1400	0,0370	3,1713	0,0373	3,2205	0,0366	3,1692	0,0406	3,2042	0,0369	3,1433	0,0389	3,2058
P	3,1008	0,0373	3,1702	0,0370	3,1252	0,0365	3,1572	0,0370	3,1393	0,0374	3,1152	0,0369	3,1270	0,0373	3,2039	0,0379	3,1550	0,0379	3,2210	0,0372	3,1727	0,0383	3,1842	0,0371

Figure 1 Absorption values at 640 nm of the microplates that were pipetted with the 50 µl SC LE (top) and 50 µl MC LE (bottom). All values greater than the critical value $A_{crit} = 0.0468$ were marked in red.

Conclusion

Pipetting with the LHS into a 384-well microplate showed no errors and no contaminations with both a single-channel LE and a multichannel LE. Therefore, the LHS is suitable for sensitive samples during HTS.



All information about our products can be found at shop.brand.de

BRAND®, BRAND. For lab. For life.® as well as the BRAND figurative mark are registered trademarks or trademarks of BRAND GMBH + CO KG, Germany. All other trademarks mentioned or depicted here are the property of the respective owners.

Our technical literature is intended to inform and advise our customers. However, the validity of general empirical values, and of results obtained under test conditions, for specific applications depends on many factors beyond our control.

Please appreciate, therefore, that no claims can be derived from our advice. The user is responsible for checking the appropriateness of the product for any particular application. California Residents: For more information concerning California Proposition 65, please refer to www.brand.de/calprop65.

Subject to technical modification without notice. Errors excepted.

BRAND GMBH + CO KG

P.O. Box 1155 | 97861 Wertheim | Germany

T +49 9342 808 0 | F +49 9342 808 98000 | info@brand.de | www.brand.de

BRAND. For lab. For life.®

