

A Closer look...

Models

Ergonomics – approved and certified.

Transferpette® electronic single channel pipette is available in 5 different models: 0.5-10 µl, 2-20 µl, 20-200 µl, 100-1000 µl and 0.5-5 ml.

The Transferpette®-8/-12 electronic multichannel pipette is available in 5 different models: 0.5-10 µl, 1-20 µl, 5-100 µl, 10-200 µl and 15-300 µl.

Designed for comfort in your relaxed hand, the Transferpette® electronic adds simple operating software and well-illustrated documentation to produce a uniquely user-friendly experience.

