According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® classic concentrate mucasolNo Change Service!VersionRevision Date:Date of last issue: 22.03.2022

Version Revision Date: 06.01 20.09.2022

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Trade name : perform® classic concentrate mucasol

1.2 Relevant identified uses of the substance or mixture and uses advised against

Use of the Sub-

stance/Mixture

: Cleaning agent

Recommended restrictions

on use

Restricted to professional users.

1.3 Details of the supplier of the safety data sheet

Producer : Schülke & Mayr GmbH

Robert-Koch-Str. 2

22851 Norderstedt

Germany

Telephone: +49 (0)40/ 52100-0 Telefax: +49 (0)40/ 52100318

mail@schuelke.com www.schuelke.com

Supplier : Schülke & Mayr UK Ltd.

Cygnet House

1, Jenkin Road, Meadowhall

Sheffield S9 1AT United Kingdom

Telephone: +44 114 254 35 00 Telefax: +44 114 254 35 01 mail.uk@schulke.com

E-mail address of person

responsible for the SDS/Contact person

Application Specialists +49 (0)40/ 521 00 666 AD@schuelke.com

(Schülke & Mayr UK Ltd.: +44-1142543500)

1.4 Emergency telephone number

Emergency telephone num-

Carechem 24 International:+44 1235 239670

ber

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Classification (REGULATION (EC) No 1272/2008) as amended by The Chemicals (Health and Safety) and Genetically Modified Organisms (Contained Use) (Amendment etc.) (EU Exit) Regulations 2019)

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® classic concentrate mucasol No Change Service! Date of last issue: 22.03.2022

06.01 20.09.2022

> Corrosive to metals, Category 1 H290: May be corrosive to metals.

Serious eye damage, Category 1 H318: Causes serious eye damage.

Specific target organ toxicity - single exposure, Category 3, Respiratory system

Revision Date:

H335: May cause respiratory irritation.

2.2 Label elements

Version

Labelling (REGULATION (EC) No 1272/2008) as amended by The Chemicals (Health and Safety) and Genetically Modified Organisms (Contained Use) (Amendment etc.) (EU Exit) Regulations 2019)

Hazard pictograms





Signal word Danger

Hazard statements H290 May be corrosive to metals.

Causes serious eye damage. H318 H335 May cause respiratory irritation.

Precautionary statements Prevention:

> P261 Avoid breathing vapours.

P280 Wear eye protection/ face protection.

Response:

P305 + P351 + P338 + P310 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a

POISON CENTER/ doctor.

Hazardous components which must be listed on the label:

tripotassium orthophosphate

2.3 Other hazards

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

SECTION 3: Composition/information on ingredients

3.2 Mixtures

Chemical nature Solution of the following substances with harmless additives.

Components

Chemical name	CAS-No.	Classification	Concentration
	EC-No.		(% w/w)
	Index-No.		

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® classic concentrate mucasolNo Change Service!VersionRevision Date:Date of last issue: 22.03.2022

Version Revision Date: 06.01 20.09.2022

Registration number		
7778-53-2 231-907-1 	Eye Dam. 1; H318 STOT SE 3; H335 (Respiratory system)	>= 20 - < 30

For explanation of abbreviations see section 16.

Other information

7778-53-2 REACH registration numbers can not be provided because the product is exempted from the registration requirements according to REACH Title II., This product is an ionic mixture.

SECTION 4: First aid measures

4.1 Description of first aid measures

General advice : Take off all contaminated clothing immediately.

If inhaled : Move to fresh air.

If symptoms persist, call a physician.

In case of skin contact : Wash with water and soap as a precaution.

In case of eye contact : In case of eye contact, remove contact lens and rinse imme-

diately with plenty of water, also under the eyelids, for at least

15 minutes.

Call a physician immediately.

If swallowed : Rinse mouth.

Do NOT induce vomiting. Obtain medical attention.

4.2 Most important symptoms and effects, both acute and delayed

Symptoms : irritant effects

Risks : Causes serious eye damage.

May cause respiratory irritation.

4.3 Indication of any immediate medical attention and special treatment needed

Treatment : For specialist advice physicians should contact the Poisons

Information Service. Treat symptomatically.

SECTION 5: Firefighting measures

5.1 Extinguishing media

Suitable extinguishing media : Carbon dioxide (CO2)

Foam Dry powder

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® classic concentrate mucasol No Change Service!

Revision Date: Date of last issue: 22.03.2022 Version

06.01 20.09.2022

Water spray jet

Unsuitable extinguishing

media

Do NOT use water jet.

5.2 Special hazards arising from the substance or mixture

Specific hazards during fire-

fighting

No information available.

ucts

Hazardous combustion prod- : No hazardous combustion products are known

5.3 Advice for firefighters

for firefighters

Special protective equipment : In the event of fire, wear self-contained breathing apparatus.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Personal precautions Increased risk of slipping in the presence of leaked / spilled

product.

Use personal protective equipment.

6.2 Environmental precautions

Environmental precautions Avoid subsoil penetration.

6.3 Methods and material for containment and cleaning up

Wipe up with absorbent material (e.g. cloth, fleece). Methods for cleaning up

Soak up with inert absorbent material (e.g. sand, silica gel,

acid binder, universal binder, sawdust).

6.4 Reference to other sections

see Section 8 + 13

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Advice on safe handling Avoid contact with skin and eyes.

Advice on protection against :

fire and explosion

No special protective measures against fire required.

Hygiene measures Keep away from food and drink.

7.2 Conditions for safe storage, including any incompatibilities

Requirements for storage areas and containers

Store at room temperature in the original container. Do not

store at temperatures above 30°C.

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® classic concentrate mucasol No Change Service! Date of last issue: 22.03.2022

06.01 20.09.2022

Further information on stor-

age conditions

Version

Keep away from heat. Keep away from direct sunlight. Keep container tightly closed. Recommended storage temperature:

5 - 25°C

Advice on common storage No materials to be especially mentioned.

7.3 Specific end use(s)

Specific use(s) none

Revision Date:

SECTION 8: Exposure controls/personal protection

8.1 Control parameters

Contains no substances with occupational exposure limit values.

Derived No Effect Level (DNEL):

Substance name	End Use	Exposure routes	Potential health effects	Value
β-Alanine, N-C8-18- alkyl derivs., monopo- tassium salts	Workers	Inhalation	Systemic effects	97.8 mg/m3
	Workers	Dermal	Systemic effects	13.9 mg/kg bw/day

8.2 Exposure controls

Personal protective equipment

Eye/face protection Safety glasses with side-shields conforming to EN166

Hand protection

Directive The selected protective gloves have to satisfy the specifica-

tions of Regulation (EU) 2016/425 and the standard EN 374

derived from it.

Splash protection: disposable nitrile rubber gloves e.g. Remarks

Dermatril (layer thickness: 0.11 mm) made by KCL or gloves from other manufacturers offering the same protection. Prolonged contact: Nitrile rubber gloves e.g. Camatril (>480 Min., layer thickness: 0,40 mm) or butyl rubber gloves e.g. Butoject (>480 Min., layer thickness: 0,70 mm) made by KCL or gloves from other manufacturers offering the same protection. Prolonged contact: Butyl rubber gloves e.g. Butoject (>480 Min., layer thickness: 0.70 mm) made by KCL or gloves from other manufacturers offering the same protec-

tion.

Skin and body protection Work uniform or laboratory coat.

Respiratory protection No personal respiratory protective equipment normally re-

quired.

Protective measures Avoid contact with eyes.

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® classic concentrate mucasol No Change Service!

Version Revision Date: Date of last issue: 22.03.2022

06.01 20.09.2022

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

Appearance : liquid

Colour : nearly colourless, -, light yellow

Odour : characteristic

Odour Threshold : not determined

pH : 12.3 (20 °C)

Concentration: 100 %

Melting point/freezing point : < 0 °C

Decomposition temperature No data available

Boiling point/boiling range : > 100 °C

Flash point : Not applicable

Evaporation rate : Not applicable

Upper explosion limit / Upper

flammability limit

No data available

Lower explosion limit / Lower

flammability limit

No data available

Vapour pressure : Not applicable

Relative vapour density : Not applicable

Density : ca. 1.43 g/cm3 (20 °C)

Solubility(ies)

Water solubility : completely miscible (20 °C)

Auto-ignition temperature : No data available

Viscosity

Viscosity, dynamic : ca. 6 mPa*s

Explosive properties : No data available

Oxidizing properties : The substance or mixture is not classified as oxidizing.

9.2 Other information

Flammability (liquids) : Does not sustain combustion.

Metal corrosion rate : > 6.25 mm/a

Corrosive to metals Aluminium and Mild steel

Z11832 ZSDB_P_GB EN

Page 6/14

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® classic concentrate mucasol No Change Service!

Version Revision Date: Date of last issue: 22.03.2022

06.01 20.09.2022

SECTION 10: Stability and reactivity

10.1 Reactivity

No dangerous reaction known under conditions of normal use.

10.2 Chemical stability

The product is chemically stable.

10.3 Possibility of hazardous reactions

Hazardous reactions : Reaction with oxidising agents

10.4 Conditions to avoid

Conditions to avoid : Protect from frost, heat and sunlight.

10.5 Incompatible materials

Materials to avoid : Possible incompatibility with alkali sensitive materials.

10.6 Hazardous decomposition products

None reasonably foreseeable.

SECTION 11: Toxicological information

11.1 Information on toxicological effects

Acute toxicity

Not classified based on available information.

Components:

tripotassium orthophosphate:

Acute oral toxicity : Remarks: No data available

Acute inhalation toxicity : Target Organs: Mucous membranes

Symptoms: Local irritation

Acute dermal toxicity : Remarks: No data available

Skin corrosion/irritation

Not classified based on available information.

Components:

tripotassium orthophosphate:

Species : reconstructed human epidermis (RhE)

Method : OECD Test Guideline 439

Result : No skin irritation

Serious eye damage/eye irritation

Causes serious eye damage.

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® classic concentrate mucasol No Change Service!

Version Revision Date: Date of last issue: 22.03.2022

06.01 20.09.2022

Components:

tripotassium orthophosphate:

Species : Rabbit

Result : Irreversible effects on the eye

Respiratory or skin sensitisation

Skin sensitisation

Not classified based on available information.

Respiratory sensitisation

Not classified based on available information.

Components:

tripotassium orthophosphate:

Test Type : Local lymph node assay (LLNA)

Species : Mouse

Method : OECD Test Guideline 429
Result : Not a skin sensitizer.

Germ cell mutagenicity

Not classified based on available information.

Components:

tripotassium orthophosphate:

Genotoxicity in vitro : Test Type: In vitro mammalian cell gene mutation test

Test system: mouse lymphoma cells

Method: QSAR Result: negative

Carcinogenicity

Not classified based on available information.

Components:

tripotassium orthophosphate:

Remarks : This information is not available.

Reproductive toxicity

Not classified based on available information.

Components:

tripotassium orthophosphate:

Effects on foetal develop- : Species: Rat

ment General Toxicity Maternal: NOAEL: > 282 mg/kg bw/day

Teratogenicity: NOAEL: > 282 mg/kg bw/day

STOT - single exposure

May cause respiratory irritation.

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® classic concentrate mucasol No Change Service!

Revision Date: Date of last issue: 22.03.2022 Version

06.01 20.09.2022

Components:

tripotassium orthophosphate:

Target Organs Respiratory system

Assessment May cause respiratory irritation.

STOT - repeated exposure

Not classified based on available information.

Aspiration toxicity

Not classified based on available information.

Further information

Product:

Remarks No data is available on the product itself.

SECTION 12: Ecological information

12.1 Toxicity

Components:

tripotassium orthophosphate:

Toxicity to fish LC50 (Oncorhynchus mykiss (rainbow trout)): > 100 mg/l

Exposure time: 96 h

Method: OECD Test Guideline 203

aquatic invertebrates

Toxicity to daphnia and other : EC50 (Daphnia magna (Water flea)): > 100 mg/l

Exposure time: 48 h

Method: OECD Test Guideline 202

Toxicity to algae/aquatic

plants

EC50 (Desmodesmus subspicatus (green algae)): > 100 mg/l

Exposure time: 72 h

Method: OECD Test Guideline 201

NOEC (Desmodesmus subspicatus (green algae)): > 100 mg/l

Exposure time: 72 h

Method: OECD Test Guideline 201

12.2 Persistence and degradability

Components:

tripotassium orthophosphate:

Biodegradability Remarks: No data available

12.3 Bioaccumulative potential

Components:

tripotassium orthophosphate:

Bioaccumulation Remarks: No data available

Z11832 ZSDB_P_GB EN Page 9/14

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019

Revision Date:



perform® classic concentrate mucasol No Change Service! Date of last issue: 22.03.2022

20.09.2022 06.01

П

Version

12.4 Mobility in soil

Components:

tripotassium orthophosphate:

Mobility Remarks: No data available

12.5 Results of PBT and vPvB assessment

Product:

Assessment : This substance/mixture contains no components considered

> to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of

0.1% or higher.

12.6 Other adverse effects

Product:

Endocrine disrupting poten-

The substance/mixture does not contain components considered to have endocrine disrupting properties according to REACH Article 57(f) or Commission Delegated regulation

(EU) 2017/2100 or Commission Regulation (EU) 2018/605 at

levels of 0.1% or higher.

Additional ecological infor-

mation

No data is available on the product itself.

SECTION 13: Disposal considerations

13.1 Waste treatment methods

Product Disposal together with normal waste is not allowed. Special

disposal required according to local regulations.

Contaminated packaging Empty containers should be taken to an approved waste han-

dling site for recycling or disposal.

SECTION 14: Transport information

14.1 UN number

ADR UN 1719 **IMDG** UN 1719 IATA UN 1719

14.2 UN proper shipping name

ADR CAUSTIC ALKALI LIQUID. N.O.S.

(Tripotassium orthophosphate)

Z11832 ZSDB_P_GB EN

Page 10/14

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® classic concentrate mucasol No Change Service!

Version Revision Date: Date of last issue: 22.03.2022

06.01 20.09.2022

IMDG : CAUSTIC ALKALI LIQUID, N.O.S.

(Tripotassium orthophosphate)

IATA : Caustic alkali liquid, n.o.s.

(Tripotassium orthophosphate)

14.3 Transport hazard class(es)

 ADR
 : 8

 IMDG
 : 8

 IATA
 : 8

14.4 Packing group

ADR

Packing group : III
Classification Code : C5
Hazard Identification Number : 80
Labels : 8
Tunnel restriction code : (E)

IMDG

Packing group : III Labels : 8

EmS Code : F-A, S-B

IATA (Cargo)

Packing instruction (cargo : 856

aircraft)

Packing instruction (LQ) : Y841
Packing group : III
Labels : Corrosive

IATA (Passenger)

Packing instruction (passen- : 852

ger aircraft)

Packing instruction (LQ) : Y841
Packing group : III

Labels : Corrosive

14.5 Environmental hazards

ADR

Environmentally hazardous : no

IMDG

Marine pollutant : no

14.6 Special precautions for user

The transport classification(s) provided herein are for informational purposes only, and solely based upon the properties of the unpackaged material as it is described within this Safety Data Sheet. Transportation classifications may vary by mode of transportation, package sizes, and variations in regional or country regulations.

14.7 Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable for product as supplied.

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® classic concentrate mucasol No Change Service!

Revision Date: Version 06.01 20.09.2022

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

Relevant EU provisions transposed through retained EU law

UK REACH List of restrictions (Annex 17) Conditions of restriction for the fol-

lowing entries should be considered:

Number on list 3

UK REACH Candidate list of substances of very high

concern (SVHC) for Authorisation

Not applicable

Not applicable

Date of last issue: 22.03.2022

The Persistent Organic Pollutants Regulations (retained Regulation (EU) 2019/1021 as amended for Great Brit-

ain)

Regulation (EC) No 1005/2009 on substances that de-

plete the ozone layer

Not applicable

UK REACH List of substances subject to authorisation Not applicable

(Annex XIV)

Directive 2010/75/EU of 24 November 2010 on industrial Volatile organic compounds

> emissions (integrated pollution prevention and control) Volatile organic compounds (VOC) content: 0.24 %

according to Detergents

30 % and more: Phosphates

Regulation EC 648/2004 less than 5 %: Amphoteric surfactants

The components of this product are reported in the following inventories:

TCSI On the inventory, or in compliance with the inventory

TSCA Product contains substance(s) not listed on TSCA inventory.

AIIC Not in compliance with the inventory

DSL This product contains the following components that are not

on the Canadian DSL nor NDSL.

β-Alanine, N-C8-18-alkyl derivs., monopotassium salts

ENCS Not in compliance with the inventory

ISHL Not in compliance with the inventory

KECI Not in compliance with the inventory

PICCS Not in compliance with the inventory

IECSC Not in compliance with the inventory

NZIoC Not in compliance with the inventory

TECI Not in compliance with the inventory

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® classic concentrate mucasol No Change Service!

Version Revision Date: Date of last issue: 22.03.2022

06.01 20.09.2022

15.2 Chemical safety assessment

Exempt

SECTION 16: Other information

Full text of H-Statements

H318 : Causes serious eye damage. H335 : May cause respiratory irritation.

Full text of other abbreviations

Eye Dam. : Serious eye damage

STOT SE : Specific target organ toxicity - single exposure

ADN - European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways; ADR - Agreement concerning the International Carriage of Dangerous Goods by Road; AIIC - Australian Inventory of Industrial Chemicals; ASTM - American Society for the Testing of Materials; bw - Body weight; CLP - Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECHA -European Chemicals Agency; EC-Number - European Community number; ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships: n.o.s. - Not Otherwise Specified: NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RID - Regulations concerning the International Carriage of Dangerous Goods by Rail; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; SVHC - Substance of very high concern; TCSI - Taiwan Chemical Substance Inventory; TECI -Thailand Existing Chemicals Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative

Further information

Classification of the mixture:

Classification procedure:

Met. Corr. 1 H290 Based on product data or assessment

Eye Dam. 1 H318 Calculation method

According to REACH etc. (Amendment etc.) (EU Exit) Regulations 2019



perform® classic concentrate mucasolNo Change Service!VersionRevision Date:Date of last issue: 22.03.2022

Version Revision Date: 06.01 20.09.2022

STOT SE 3 H335 Calculation method

Changes since the last version are highlighted in the margin. This version replaces all previous versions.

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.