



Microplate, 384-well 781680



L	127.55 mm
W	85.05 mm
H	11.50 mm
Hs	3.50 mm
X1	12.03 mm
Xn	115.53 mm
Y1	8.78 mm
Yn	76.28 mm
P	4.50 mm
D	3.90 mm
Hw	8.85 mm
C	24 Well
R	16 Well




Type	standard
Base shape	F-bottom
Surface	pureGrade™ S
Volume	100 µl / well
Material	PS
Color	transparent
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, Sterile according to ISO 11 137, Free of cytotoxic substances acc. ISO 10 993-5




* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml

Subject to technical modification without notice. Errors excepted.

Microplate, 384-well 781681



L	127.55 mm
W	85.05 mm
H	11.50 mm
Hs	3.50 mm
X1	12.03 mm
Xn	115.53 mm
Y1	8.78 mm
Yn	76.28 mm
P	4.50 mm
D	3.90 mm
Hw	8.85 mm
C	24 Well
R	16 Well




Type	standard
Base shape	F-bottom
Surface	pureGrade™ S
Volume	100 µl / well
Material	PS
Color	white
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, Sterile according to ISO 11 137, Free of cytotoxic substances acc. ISO 10 993-5




* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml

Subject to technical modification without notice. Errors excepted.

Microplate, 384-well 781682



L	127.55 mm
W	85.05 mm
H	11.50 mm
Hs	3.50 mm
X1	12.03 mm
Xn	115.53 mm
Y1	8.78 mm
Yn	76.28 mm
P	4.50 mm
D	3.90 mm
Hw	8.85 mm
C	24 Well
R	16 Well



Type	standard
Base shape	F-bottom
Surface	pureGrade™ S
Volume	100 µl / well
Material	PS
Color	black
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, Sterile according to ISO 11 137, Free of cytotoxic substances acc. ISO 10 993-5




* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml

Subject to technical modification without notice. Errors excepted.

Microplate, 384-well

781686



L	127.60 mm
W	85.60 mm
H	14.80 mm
Hs	2.75 mm
X1	12.10 mm
Xn	115.60 mm
Y1	9.10 mm
Yn	76.60 mm
P	4.50 mm
D	3.65 mm
Hw	11.95 mm
C	24 Well
R	16 Well


Type	standard, transparent bottom
Base shape	F-bottom
Surface	pureGrade™ S
Volume	120 µl / well
Material	PS
Color	white
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, Sterile according to ISO 11 137, Free of cytotoxic substances acc. ISO 10 993-5



* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml

Subject to technical modification without notice. Errors excepted.

Microplate, 384-well 781687



L	127.60 mm
W	85.60 mm
H	14.80 mm
Hs	2.75 mm
X1	12.10 mm
Xn	115.60 mm
Y1	9.10 mm
Yn	76.60 mm
P	4.50 mm
D	3.65 mm
Hw	11.95 mm
C	24 Well
R	16 Well

Type	standard, transparent bottom
Base shape	F-bottom
Surface	pureGrade™ S
Volume	120 µl / well
Material	PS
Color	black
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, Sterile according to ISO 11 137, Free of cytotoxic substances acc. ISO 10 993-5



* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml

Subject to technical modification without notice. Errors excepted.