


# Single Insert, height 13 mm for BRANDplates® 24-well plate 782806



PRODUCT  
INFORMATION



Dimensions in mm

<b>Volume</b>	50 - 400 µl
<b>Material</b>	Polystyrene, medical grade
<b>Quality</b>	BIO-CERT® CELL CULTURE QUALITY




## Membrane properties

<b>Material</b>	Polycarbonate, USP class VI conform
<b>Thickness</b>	10 µm
<b>Pore size</b>	0.4 µm
<b>Pore density</b>	1 x 10 <sup>8</sup> /cm <sup>2</sup> , translucent
<b>Culture area</b>	ca. 0.6 cm <sup>2</sup>
<b>Surface treatment</b>	cellGrade™ plus

Subject to technical modification without notice. Errors excepted.

BRAND GMBH + CO KG | Otto-Schott-Str. 25 | 97877 Wertheim | Germany | T +49 9342 808 0 | F +49 9342 808 98000 | info@brand.de | www.brand.de

# Single Insert, height 13 mm for BRANDplates® 24-well plate 782826


 PRODUCT  
INFORMATION


Dimensions in mm

<b>Volume</b>	50 - 400 µl
<b>Material</b>	Polystyrene, medical grade
<b>Quality</b>	BIO-CERT® CELL CULTURE QUALITY




## Membrane properties

<b>Material</b>	Polycarbonate
<b>Thickness</b>	9 µm
<b>Pore size</b>	1.0 µm
<b>Pore density</b>	2x 10 <sup>6</sup> /cm <sup>2</sup>
<b>Culture area</b>	ca. 0.6 cm <sup>2</sup>
<b>Surface treatment</b>	cellGrade™ plus

Subject to technical modification without notice. Errors excepted.

BRAND GMBH + CO KG | Otto-Schott-Str. 25 | 97877 Wertheim | Germany | T +49 9342 808 0 | F +49 9342 808 98000 | info@brand.de | www.brand.de

# Single Insert, height 13 mm for BRANDplates® 24-well plate 782846


 PRODUCT  
INFORMATION


Dimensions in mm

<b>Volume</b>	50 - 400 µl
<b>Material</b>	Polystyrene, medical grade
<b>Quality</b>	BIO-CERT® CELL CULTURE QUALITY




## Membrane properties

<b>Material</b>	Polycarbonate
<b>Thickness</b>	22 µm
<b>Pore size</b>	3.0 µm
<b>Pore density</b>	1,6 x 10 <sup>6</sup> /cm <sup>2</sup>
<b>Culture area</b>	ca. 0.6 cm <sup>2</sup>
<b>Surface treatment</b>	cellGrade™ plus

Subject to technical modification without notice. Errors excepted.

BRAND GMBH + CO KG | Otto-Schott-Str. 25 | 97877 Wertheim | Germany | T +49 9342 808 0 | F +49 9342 808 98000 | info@brand.de | www.brand.de

# Single Insert, height 13 mm for BRANDplates® 24-well plate 782866


 PRODUCT  
INFORMATION


Dimensions in mm

<b>Volume</b>	50 - 400 µl
<b>Material</b>	Polystyrene, medical grade
<b>Quality</b>	BIO-CERT® CELL CULTURE QUALITY



## Membrane properties

<b>Material</b>	Polycarbonate
<b>Thickness</b>	18 µm
<b>Pore size</b>	8.0 µm
<b>Pore density</b>	1 x 10 <sup>5</sup> /cm <sup>2</sup>
<b>Culture area</b>	ca. 0.6 cm <sup>2</sup>
<b>Surface treatment</b>	cellGrade™ plus

Subject to technical modification without notice. Errors excepted.

BRAND GMBH + CO KG | Otto-Schott-Str. 25 | 97877 Wertheim | Germany | T +49 9342 808 0 | F +49 9342 808 98000 | info@brand.de | www.brand.de