


# ASTM Centrifuge Tubes, 100 ml

Type: cylindrical, conical bottom, length max. 203 mm

Range	Graduation [ml]	Subdivision [ml]	Error limit [± ml]
<b>BLAUBRAND®</b>			
0 - 0.1		0.05	0.02
> 0.1 - 0.3		0.05	0.03
> 0.3 - 0.5		0.05	0.05
> 0.5 - 1		0.1	0.05
> 1 - 2		0.1	0.1
> 2 - 3		0.2	0.1
> 3 - 5		0.5	0.2
> 5 - 10		1	0.5
> 10 - 25		5	1
> 25 - 100		25	1



Type: cylindrical, conical bottom, length max. 167 mm

Range	Graduation [ml]	Subdivision [ml]	Error limit [± ml]
<b>BLAUBRAND®</b>			
0 - 0.1		0.05	0.02
> 0.1 - 0.3		0.05	0.03
> 0.3 - 0.5		0.05	0.05
> 0.5 - 1		0.1	0.07
> 1 - 1.5		0.1	0.1
> 1.5 - 2		0.1	0.2
> 2 - 3		0.2	0.3
> 3 - 5		0.5	0.5
> 5 - 10		1	0.75
> 10 - 25		5	1
> 25 - 100		at 50 and 100	1.5

Type: pear-shaped, cylindrical bottom, length max. 160 mm

Range	Graduation [ml]	Subdivision [ml]	Error limit [± ml]
<b>BLAUBRAND®</b>			
0 - 1.5		0.1	0.03
> 1.5 - 3		0.5	0.2
> 3 - 5		0.5	0.3
> 5 - 10		1	0.5
> 10 - 25		5	1
> 25 - 100		25	2