

Volumetric Flasks

Standard type

Capacity [ml]	Neck size	Sort of neck	Height [mm]	Neck inner diameter [mm]	Error limit [± ml]		
					BLAUBRAND®	BLAUBRAND® USP	SILBERBRAND
5	NS 10/19	Narrow Neck	70	7 ±1	–	0.02	–
5	NS 10/19	Wide Neck	70	9 ±1	0.04	–	–
10	NS 7/16	Narrow Neck	90	7 ±1	0.025	–	–
10	NS 10/19	Narrow Neck	90	7 ±1	–	0.02	–
10	NS 10/19	Wide Neck	90	9 ±1	0.04	–	0.06
20	NS 10/19	Narrow Neck	110	7 ±1	–	0.02	–
20	NS 10/19	Narrow Neck	110	9 ±1	0.04	–	–
20	NS 12/21	Wide Neck	105	11 ±1	0.06	–	–
25	NS 10/19	Narrow Neck	110	9 ±1	0.04	0.03	0.06
25	NS 12/21	Wide Neck	110	11 ±1	0.06	–	–
50	NS 12/21	Narrow Neck	140	11 ±1	0.06	0.05	0.09
50	NS 14/23	Wide Neck	140	13 ±1	0.10	–	–
100	NS 12/21	Narrow Neck	170	13 ±1	0.10	–	–
100	NS 14/23	Narrow Neck	170	13 ±1	0.10	0.08	0.15
150	NS 14/23	Narrow Neck	185	15.5 ±1.5	0.15	–	–
200	NS 14/23	Narrow Neck	210	15.5 ±1.5	0.15	0.10	0.25
250	NS 14/23	Narrow Neck	220	15.5 ±1.5	0.15	0.12	0.25
500	NS 19/26	Narrow Neck	260	19 ±2	0.25	0.20	0.40
1000	NS 24/29	Narrow Neck	300	23 ±2	0.40	0.30	0.60
1000	NS 29/32	Wide Neck	300	27.5 ±2.5	0.60	–	–
2000	NS 29/32	Narrow Neck	370	27.5 ±2.5	0.60	0.50	0.90
5000	NS 34/35	Narrow Neck	475	38 ±3	1.2	–	1.8
10000	NS 45/40	Narrow Neck	570	48 ±4	2.0	–	–

Trapezoidal

Capacity [ml]	Neck size	Sort of neck	Height [mm]	Neck inner diameter [mm]	Error limit [± ml]
					BLAUBRAND®
1	NS 7/16	Narrow Neck	65	7 ±1	0.025
2	NS 7/16	Narrow Neck	70	7 ±1	0.025
5	NS 7/16	Narrow Neck	70	7 ±1	0.025
10	NS 10/19	Wide Neck	90	9 ±1	0.04
20	NS 10/19	Narrow Neck	110	9 ±1	0.04
25	NS 10/19	Narrow Neck	110	9 ±1	0.04
50	NS 12/21	Narrow Neck	140	11 ±1	0.06

