

Plastic Cylinder

Graduated cylinder: tall form

Capacity [ml]	Subdivision [ml]	Height [mm]	Outer Ø [mm]	Width across corners x width across flats [mm]	Error limit [± ml]	
					class A	class B
10	0.2	145	14	50 x 43	0.10	0.20
25	0.5	170	20	55 x 48	0.25	0.5
50	1	200	26	60 x 52	0.5	1.0
100	1	250	31	70 x 61	0.5	1.0
250	2	315	42	90 x 78	1.0	2.0
500	5	360	54	110 x 95	2.5	5
1000	10	440	65	130 x 113	5	10
2000	20	535	84	159 x 139	10	20

Graduated cylinder: low form

Capacity [ml]	Subdivision [ml]	Height [mm]	Outer Ø [mm]	Diameter base [mm]	Error limit [± ml]	
					class B	
25	0.5	122	26	45	0.5	
50	1	142	30	49	1.0	
100	2	163	39	59	2.0	
250	5	192	54	74	5	
500	10	218	69	99	10	
1000	20	285	81	121	20	

